

INGEZIM PPA CROM Anticuerpo

R.11.PPA.K41

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

INGEZIM PPA CROM Anticuerpo

R.11.PPA.K41

INGEZIM PPA CROM is based on a Direct Immunochromatography technique which uses protein VP72 of ASFV.

TECHNICAL BASIS

The device consists of a plastic casing with two windows:

- Window for sample addition: Contains VP72 and a control protein, bound to coloured latex.
- Window for reading of results: Contains a test line (T) made of VP72 and a control line (C) made of an antibody specific to the control protein.

When the sample is added, if it contains specific antibodies, they will bind to the VP72 conjugated with latex and will flow along the membrane. The antigen-antibody-latex complex will bind to the antigen placed on the test zone (T) showing a pink line. The blue line in the control zone indicates the validation of the assay.

APPLICATION

Detection of specific antibodies to ASFV's VP72 in serum, plasma and full blood samples in domestic pig and wild boars.

DIAGNOSTIC SENSITIVITY

- 84 warthog sera from East Africa positive by OIE ELISA.
- 15 Spanish domestic pig positive sera by OIE ELISA.
- 109 Western African domestic pig positive sera by OIE ELISA.
- 13 Italian (Sardinia) wild boars sera positive by IPMA

99% correspondence with the OIE ELISA. 82% sensitivity respect to IPMA (wild boars)

DIAGNOSTIC SPECIFICITY

- 1043 domestic pig sera from ASFV free zones.
- 167 African domestic pig sera negative by OIE ELISA and INGEZIM® PPA COMPAC.
- 97 Italian (Sardinia) wild boars sera negative by IPMA.

99.9% correspondence with INGEZIM® PPA COMPAC and OIE ELISA. 96% specificity respect IPMA (wild boars)

ANALYTICAL SENSITIVITY

- Sera from 13 ASFV experimentally infected pigs.
- OIE Reference Sera : Strong positive and weak positive

The assay is able to detect antibodies in the range of day 10p.i. and day 21p.i. depending on the viral isolate.

Detects OIE Strong positive serum at dilution 1/64 in negative serum and the OIE Weak positive serum at dilution 1/2 in negative serum.

These studies have been made in collaboration with CISA-INIA (ASFRISK Project)

Cappai et al (2017) Evaluation of a Commercial Field Test to Detect African Swine Fever. J. of Wildlife diseases, 53(3)

COMPOSITION OF THE KIT

- Chromatographic devices
- Vials with diluent

PRODUCT MANUFACTURED BY INGENASA
REGISTRATION NUMBER 1719 RD

SHELF LIFE: **24 months**
Stored at 4°C-25°C

Ed.020217