INgezim® CIRCO IgG

R.11.PCV.K1

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
REFERENCIAS: (A) Segalés et al. 6th International Symposium on Emerging and Re-emerging Pig Diseases Barcelona, Spain, 2011); (B) Pileri et al. 22nd International Pig Veterinary Society Congress. Jeju, Korea (2012); (B)Segalés et al. AVEDILA 2013, Madrid Spain.

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INgezim® Circo IgG

R.11.PCV.K1

INgezim® CIRCO IgG is an immunoenzymatic assay based on indirect ELISA technique, which uses a monoclonal antibody (Mab) specific to porcine immunoglobulins, and a recombinant antigen (PCV2 VP2 protein).


KIT FEATURES

APPLICATION

Detection and/or titration of specific antibodies to PCV in porcine sera samples.

TECHNICAL BASE

- Plates are supplied coated with PCV2 antigen (VP2 recombinant protein). Samples are added to the wells and incubated.
- If the samples contain specific antibodies to PCV, they will bind to the antigen.
- When the conjugate (monoclonal antibody specific to porcine immunoglobulins, conjugated with peroxidase, AcM-PO) is added, it will bind to the antibodies bound to the antigen. This binding is revealed by colorimetric reaction after addition of the substrate.

RESULTS INTERPRETATION


The assay establishes a cut off, which will classify the samples as *Positive* or *Negative*, depending on the value of the optical density of the sample. The titer of the sample is the last dilution showing an OD higher than the cut off.

One well titration: it is possible to obtain an approximated titre of the sera based on its OD value, by using a formula which corresponds with an exponential curve.

ASSAY VALIDATION


Experimental study 1

A set of sera from three experimentally infected piglets previously vaccinated were used. The samples consisted of different extractions made at days 0, 16 and 46 post vaccination and 14 and 22 post infection. Piglets 40 and 70 showed presence of maternal antibodies at day 0.

Studies with piglets

Studies with piglets from vaccinated sows were performed. Extractions between 3 and 24 weeks old were taken. The results were compared with the ones obtained using INgezim® Circovirus IgG/IgM.


Correlation with IPMA

Study A: This study was made with 645 sera from animals of different ages.

Study B This study was made with 1248 sera from animals of different ages.

Log₂ PCV2 IPMA titers correspondence with ELISA OD (A) or S/P (B) values.

Lineal c. ($R^2=0,88$); Polynomial c. ($R^2=0,99$) Lineal c. ($R^2=0,71$); Polynomial c. ($R^2=0,91$)

Experimental study 2

Sera from 12 experimentally infected animals and 8 non infected animals were analysed. Samples consisted of extractions taken at days 14 and 21 post infection. The results indicated that the assay is able to detect antibodies at least at day 21 post infection.

REFERENCES: (A) Segalés et al. 6th International Symposium on Emerging and Re-emerging Pig Diseases Barcelona, Spain, 2011); (B) Pileri et al. 22nd International Pig Veterinary Society Congress. Jeju, Korea (2012); (B)Segalés et al. AVEDILA 2013, Madrid Spain.

KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

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Eurofins-INGENASA, S.A.