## INgezim® ADV TOTAL

R.11.ADV.K1

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# INgezim® ADV TOTAL

R.11.ADV.K1

INgezim® ADV TOTAL is an immunoenzymatic assay based on indirect ELISA technique, which uses a monoclonal antibody (Mab) specific to porcine IgG and inactivated ADV antigen.

## KIT FEATURES

### APPLICATION


Detection and/or titration of specific antibodies to ADV in porcine sera samples.

### TECHNICAL BASE

1. Plates are supplied coated with antigen (inactivated ADV antigen). Samples are added to the wells and incubated.
2. If the samples contain specific antibodies to ADV, they will bind to the antigen.
3. When the conjugate (monoclonal antibody specific to porcine IgG, conjugated with peroxidase, AcM-PO) is added, it will bind to the antibodies bound to the antigen. This binding is revealed by colorimetric reaction after addition of the substrate.

### RESULTS INTERPRETATION

The assay establishes two cut offs, which will classify the samples as *Positive* or *Negative*, depending on the value of the optical density of the sample in the assay.

## ASSAY VALIDATION

### SENSITIVITY AND SPECIFICITY

#### 1. Using reference sera from O.I.E.

The international reference serum ADV-GI (1/8) is detected as required by the O.I.E. A set of reference sera (Q, B, G, D, J and F) catalogued as positive or negative by a O.I.E. Reference Laboratory (Maissons Alfort, France) were analysed. The correspondence with the expected results was 100%.

#### 2. Correspondence with the O.I.E. reference technique SN (seroneutralization)

- o Using field sera. 142 field sera were analysed in order to determine the sensitivity and specificity of the INgezim® ADV TOTAL assay in relation to the SN technique.
- o Using experimentally vaccinated sera. 210 samples from 18 experimentally vaccinated animals (with a deleted vaccine) and 8 negative controls were analysed. The samples were taken at different days starting at day 0, for 4 months. The results obtained indicated that the correspondence between INgezim® ADV TOTAL and Seroneutralization was 97%.

#### 3. Sensitivity using experimentally vaccinated sera

Experimentally vaccinated animals were used, and the antibodies detection profile was obtained. The assay can detect antibodies specific to ADV from day 13 post vaccination.

#### 4. Specificity using experimentally vaccinated animals with a different pathology (PRRS)

96 sera from animals experimentally vaccinated with a vaccine not related to ADV were analysed. The results obtained indicated that the presence of different antibodies does not interfere with the assay, showing a specificity of 100%.

#### 5. External validation

The product has been included in the study "Essai Interlaboratoires D' Aptitude Maladie D' Aujeszky 2003", organized by the O.I.E. Reference Centre for ADV. INgezim® ADV TOTAL obtained a correspondence of 100% with the expected results.

## KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

Spanish registration nº 392RD

EXPIRATION: 18 MONTHS. Stored at 2°C-8°C

**Eurofins-INGENASA**  
 Avda. de la Institución Libre de Enseñanza 39, 8º  
 28037 MADRID (SPAIN)  
 Phone: (+34)91 3680501  
[www.ingenasa.com](http://www.ingenasa.com)

IT-73840  
IT-73780

9191.INGE

9175.ING2