INgezim® ADV gE PLUS

R.11.GEP.K3

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INgezim® ADV gE PLUS

R.11.GEP.K3

INgezim® ADV gE PLUS is an immunoenzymatic assay based on a blocking ELISA assay technique, which uses a monoclonal antibody (Mab) specific to Aujeszky Disease Virus (ADV) gE protein.


KIT FEATURES

APPLICATION

Detection of specific antibodies to ADV gE protein in porcine sera samples. The assay allows the specific differentiation between antibodies produced by infection or vaccination. The assay has been developed to avoid specificity problems due to unspecific antibodies present in sera from breeding animals, without losing sensitivity in the test.

TECHNICAL BASE

- Plates are supplied coated with inactivated ADV antigen. Serum samples are added to the wells and incubated.
- If the samples contain specific antibodies to ADV, they will bind to the antigen.
- When the conjugate (monoclonal antibody specific to ADV gE protein, conjugated with peroxidase, AcM-PO) is added, it will bind to the protein only if there are no antibodies in the sample blocking the antigen (negative animals). In case there were antibodies blocking the antigen (infected animals), the conjugate would not be able to bind it. This binding is revealed by colorimetric reaction after addition of the substrate.

RESULTS INTERPRETATION

The assay establishes two cut offs, which will classify the samples as **Positive** or **Negative**, depending on the value of the optical density of the sample in the assay.

ASSAY VALIDATION

SENSITIVITY AND SPECIFICITY

- Using reference sera from O.I.E.**
The international reference serum ADV-GI (1/8) and the equivalent one are detected as required by the O.I.E.
- Specificity using sera from multi-vaccinated animals and from animals infected with different pathogens other than ADV.**
A set of 495 sera from multi-vaccinated animals with different deleted vaccines (GESKIPUR, NOVIPORVAC, SUVAXYN, AUGIVAC GP, NIA 4 and BEGONIA), and 40 sera from experimentally infected animals (with ASFV, Erysipelas, TGEV and PRRS) were analysed with INgezim® ADV gE. The results obtained were negative in 100% of the sera samples.
- Specificity using sera from breeding animals.**
509 sera from breeding animals obtained from 11 different problematic herds were analysed. The results obtained showed **99.8% specificity** and **98.9% sensitivity** when these kinds of sera are analysed.
- Specificity using sera from experimentally vaccinated animals.**
Sera of 18 experimentally vaccinated animals with a deleted vaccine were analysed. In the same study, sera of 8 non vaccinated animals were checked as controls. The results obtained indicated that the assay was 100% specific as no positive results were obtained.
- External validation**
This product has been included in the "SEROAUJ 2015", organized by the OIE Reference Laboratory (ANSES, France). In this report it is described that INgezim® ADV gE PLUS is able to detect a 1/128 dilution of the serum S96AUJ02 which indicates that it is able to detect also as positive the dilution 1/64 which is equivalent to dilution 1/8 of the OIE Reference Serum ADV-1.
INgezim® ADV gE PLUS showed 100% correspondence with the expected results in this study.

KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

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