## INgezim® FVR Compac

R.13.FVR.K3

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# INgezim® FVR Compac

R.13.FVR.K3

INgezim® FVR Compac is an immunoenzymatic assay based on a blocking ELISA technique, which uses a monoclonal antibody (MAb) specific to Rift Valley Fever Virus (RVFV) N protein and the recombinant N protein of the RVFV.


## KIT FEATURES

### APPLICATION

Detection of specific antibodies to the N protein of RVFV. The detection of antibodies does not depend on the animal species and can therefore be used with sera samples of all species affected by RVFV.

### TECHNICAL BASE

- Plates are supplied coated with antigen (recombinant N protein of RVFV). Samples are added to the wells and incubated.
- If the samples contain specific antibodies to RVFV, they will bind to the antigen.
- When the conjugate (monoclonal antibody specific to RVFV N protein, conjugated with peroxidase, AcM-PO) is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the antigen. In case the sample contains antibodies blocking the antigen (infected animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

### RESULTS INTERPRETATION

The assay establishes two cut off, which will classify the samples as *Positive* or *Negative*, depending on the Inhibition Percentage (IP) value of the sample in the assay, considering a range of inhibition percentages close to the cut off as *Doubtful* results.

## ASSAY VALIDATION

### SENSITIVITY

The following were analysed:


- 23 sheep experimentally infected with RVFV and bled at days 18 and 24 post infection.
- 8 sheep infected with RVFV.

The results obtained indicated 97% sensitivity.

### SPECIFICITY

1526 animals (sheep, goats, and cattle) from RVFV free areas were analysed. The results obtained indicated specificity values higher than 99%.

On the other hand, 1014 wild animals (deer, ibex, mouflons, fallow deer, alpaca, and zebra) from RVFV free areas in Sapin were analysed, being the specificity of the assay higher than 99%.

## KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (ABTS)
- Bottles with Stop Solution

Spanish registration nº 3300RD

EXPIRATION: 15 MONTHS. Stored at 2°C-8°C

Eurofins-INGENASA, S.A

Avda. de la Institución Libre de Enseñanza 39, 8º  
28037 MADRID (SPAIN)  
Phone: (+34)91 3680501  
[www.ingenasa.com](http://www.ingenasa.com)

IT-73840  
IT-73780

9191.INGE

9175.ING2