INgezim® ADV gB

R.11.GB.K3

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INgezim® ADV gB

R.11.GB.K3

INgezim® ADV gB is an enzymatic assay based on a blocking ELISA technique, which uses a monoclonal antibody (MAb) specific to Aujeszky Disease Virus (ADV) gB protein.


KIT FEATURES

APPLICATION

Detection of specific antibodies to Aujeszky disease virus gB glycoprotein in pig serum.

TECHNICAL BASE

- Plates are supplied coated with inactivated ADV antigen. Samples are added to the wells and incubated.
- If the samples contain specific antibodies to ADV, they will bind to the antigen.
- When the conjugate (monoclonal antibody specific to ADV gB protein, conjugated with peroxidase, AcM-PO) is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the protein. In case the sample contains antibodies blocking the antigen (infected or vaccinated animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

RESULTS INTERPRETATION

The assay establishes one cut off, which will classify the samples as **Positive** or **Negative**, depending on the value of the inhibition percentage of the sample in the assay.

ASSAY VALIDATION

Specificity of the assay in breeding farms

663 swine sera from 22 breeding farms with A4 status regarding to ADV were analysed.

The results that were obtained indicated 99.7% specificity of the assay, when negative samples from breeding farms are analysed.

Sensitivity of the assay in vaccinated animal farms

Sera from 362 vaccinated animals from 21 farms were analysed, classified with respect to their status in relation to ADV as A3, both from fattening and breeding animals.

The results that were obtained indicated 99.6 sensitivity, when samples from vaccinated animals are analysed.

Sensitivity and specificity of the assay in non-breeding farms


353 sera from non-breeding sows were analysed, including sera from vaccinated fattening pigs from 4 farms and sera from ADV-negative pigs from 16 farms, both from fattening animals and boars.

The results that were obtained in this study showed specificity and sensitivity greater than 99.9%, when sera from non-breeding farms are analysed.

Global performance of the assay

A total of 1300 samples from 59 farms were analysed, including samples of breeding sows from A4 and A3 farms, from vaccinated and non-vaccinated fattening animals, as well as samples from boars from A4 farms.

The results that were obtained in this study showed 99.8% specificity and 99.6% sensitivity.

KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

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Eurofins-INGENASA, S.A