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
I. TECHNICAL BASIS

The gB test uses a monoclonal antibody (Mab) specific for gB glycoprotein of ADV and microtiter plates coated with ADV antigen. After sample addition, if antibodies are absent in the test serum, the gB antigen will remain free and the Mab peroxidase-conjugated specific for gB, will bind. After substrate addition a colour

reaction is produced. On the other hand, if gB antibodies are present in the test sample they will bind to the gB antigen of the well and will block the binding of the conjugated Mab. In this case no colour reaction will be produced.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handle
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.
10. Substrate must be handle with care, it is very sensible to light and contamination
11. Stop solution is a strong acid. Handle with care.

III. STORAGE OF COMPONENTS

All reagents and components must be stored between +2°C and +8°C. **Once opened**, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

Warm a 37°C the serum diluent if precipitation occurs.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLE

Serum samples and controls must be used without previous dilution adding 80µl/well.

VI. PREPARATION OF REAGENTS

• **Washing solution:**

Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionized water. When ready this solution remains stable when is stored between +2°C and +8°C.

• **Conjugate:**

The conjugate is ready to use.

♦ **Controls (+) & (-) :**

Control sera should be used like samples (80 µl/well).

VII. TEST PROCEDURE

- All reagents must be allowed to come to room temperature before use.
- Add 20 µl of serum diluent to each well and then, 80 µl of serum (controls and samples) to each well. Gently, mix the plate, seal it and incubate for **1 hour at 37°C or 16-18 hours at 4°C**.
- Wash 3 times following the described procedure (point IV)
- Add 100 µl of conjugate to each well. Seal the plate and **incubate for 30 minutes at 37°C**
- Wash 5 times following the procedure.
- Add 100 µl of substrate, to each well. Keep the plate for **10 min at room temperature** in a dark place.
- Add 100 µl of stop solution to each well.
- Read the OD at 450 nm within 5 min after the addition of stop solution

VIII. READING AND RESULT INTERPRETATION

A.- TEST VALIDATION

The test must be considered valid if:

- **Positive control OD is Lower than 0.3**
- **Negative control OD is higher than 0.8**

B.- INHIBITION PERCENTAGE (IP) VALUES CALCULATION:

$$IP = [1 - (M / CN)] \times 100$$

Where: CN= Mean of OD_{450nm} of Negative control
M=Mean of OD_{450nm} of the sample

C.- RESULT INTERPRETATION:

IP higher than 30%: Positive samples
IP lower than 30%: Negative samples

D.- PRACTICAL EXAMPLE:

OD₄₅₀ CN = 1.360
OD₄₅₀ CP = 0.110
OD₄₅₀ Sample 1 = 0.373
OD₄₅₀ Sample 2 = 1.295

Validation:
CN = 1.360 (> 0.8)
CP = 0,110 (< 0.3)

IP sample 1 = [1-(0.373/1.360)]x100=72.5
IP sample 2 = [1-(1.295/1.360)]x100=4.7

Interpretation:

Sample 1: Positive to gB antibodies (the animal has been in touch with a vaccinal and/or field ADV strain).

Sample 2: Negative to gB antibodies (the animal has never been in touch with any of ADV strains. This result will be obtained with non-infected and non-vaccinated pigs).

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

Eurofins-INGENASA
Av.de la Institución Libre de Enseñanza 39,
8ª planta
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in by: