INgezim® ASFV CROM Ag

R.11.ASF.K.42

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INgezim® ASFV CROM Ag

R.11.ASF.K.42

INgezim® ASFV CROM Ag is an enzymatic assay based on a direct immunochromatographic technique in which a monoclonal antibody (MAb) specific of VP72 of African Swine Fever virus (ASFV) conjugated with latex microspheres is used.

KIT FEATURES


APPLICATION

Detection of African Swine Fever Virus (antigen) in blood samples of domestic pigs and wild boars.

TECHNICAL BASE

This diagnostic device consists on a plastic plate containing a membrane. On the membrane, two lines are printed: the test line (T) has a specific MAb to ASFV and the control line (C) has a specific MAb for the control protein.

If the case of a positive sample, the virus binds to the black beads conjugated to the anti-ASFV MAb. The immune complex then migrates through the membrane by capillarity and is captured again by the anti-ASFV MAb absorbed on the test line, resulting in the appearance of a black line. The presence of the control line serves as validity of the test, indicating that the immunochromatography has been performed correctly.

ASSAY VALIDATION

DIAGNOSTIC SENSITIVITY AND SPECIFICITY

Study 1. This assay was evaluated with field and experimental samples all of them previously classified by rt-PCR Universal Probe Library (UPL) as "Gold Standard" :

- 153 blood samples (60 positive & 93 negative) of 30 pigs. These pigs were experimentally infected with different genotypes (n=23, IX; n=19, Ucrania II y n=111, Lituania II)* and bled at different days post infection.
- 58 field samples (52 positive y 6 negative) of domestic pigs and wild boards from ASFV affected countries during 2014 and 2015.*
- 15 blood samples (12 positive y 3 negative) of domestic pigs and wild boars from areas affected with genotype II.

36 samples gave false negative results. 20 of these samples gave Ct values higher than 30 and were samples taken off at early states of infection or corresponded to animals in contact with LY40/1490 isolate exposed animals. 10 of these 20 samples showed negative results by isolation. Comparing the results with the obtained by INgezim® PPA DAS (Ag detection ELISA) the coefficient κ was 0.92, indicating a perfect coincidence between both techniques.

Study 2. The assay was evaluated with 42 (37 positive y 5 negative) field samples at Friedrich Loeffler Institute (FLI), all of them previously characterized by rt-PCR. INgezim® ASF CROM Ag detected 39 positive samples and 3 negative ones.

From both studies it can be concluded that the sensitivity of the assay respect to rt-PCR was 76% and the specificity 96%.


ANALYTICAL SENSITIVITY

To determine the analytical sensitivity, blood negative samples contaminated with recombinant vp72 or with inactivated virus (strain BA71) were used. Serial dilutions of both kind of matrixes were made and the ability of detection were tested in comparison with INgezim® PPA DAS.

INgezim® ASF CROM Ag was able to detect 3ng/test of recombinant purified VP72 and less than 104 ufp/ml of viral isolate.

All these studies have been made within the frame of the project ASFORCE in collaboration with CISA-INIA.

*Sastre et al (2016) Development of a novel LFA for detection of African swine fever in blood. *Vet. Research*, 12:206

COMPOSITION OF THE KIT

- Chromatographic devices
- Vials with diluent

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EUROFINS-INGENASA, S.A.