INgezim PRRS 2.0

Prod Ref: 11.PR2.K1

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Enzymatic immunoassay for the detection and quantification of antibodies specific for porcine respiratory and reproductive syndrome virus (European and American strains) in swine serum

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Changes on the insert have been made.

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Distilled or deionised water.
Micropipets from 5 to 200 µl.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 ml
ELISA Reader (450 nm filter)

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I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). The technique is briefly described below:

The recombinant antigen, specific for antibodies of both European and American strains of PRRS virus, is fixed on polystyrene plates. When a serum sample contains specific antibodies against the virus, they will bind to the recombinant antigen adsorbed on the plate. After a washing step to eliminate all non-fixed material, the presence of swine antibodies can be detected using a specific peroxidase

conjugate. After the addition of the substrate a colorimetric reaction will appear which can be measured by a spectrophotometer at 450nm.

Please note that our kit uses purified antigen obtained by the expression of the ORF7 from American and European strains of PRRS virus. This method warrants the total absence of infectivity in any of the components of the kit.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions for use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix reagents or use instructions from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiry dates and do not mix components from different batches.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Systematically include a positive control and a negative serum each time the assay is run.
10. Substrate must be handled with care, as it is very sensitive to light and contamination.
11. Stop solution is a strong acid and therefore must be handled with care.

III. STORAGE OF COMPONENTS

Store all plates and reagents between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING PROCEDURE

The washing steps can be done using an automatic plate washer or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the instructions below:

- Throw away the content of the plate by a briskly turning the plate over in order to avoid the possible exchange of the contents from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.

- Shake the plate gently, avoiding cross-contamination between wells.
- Briskly turn the plate over to empty the wells.
- Repeat the process as many times as indicated in the kit's instructions.
- Prior to emptying the content of the plate in the last washing step, verify that the next reagent to be added is ready to used. Do not let the plate dry longer than strictly necessary.
- After the last washing step tap the plate upside down on absorbent filter paper to remove any remaining washing solution.

V. PREPARATION OF SAMPLES

Sera samples μ must be tested at **1/40 dilution in diluent DE35-01**

(i.e.: 5 μ l of serum + 195 μ l of serum diluent DE35-01).

VI. PREPARATION OF REAGENTS

- **Washing solution:**
Dilute one part of the concentrated washing solution provided in the kit with 24 parts of distilled or deionized water (40 ml of concentrated solution and 960ml of water). Once prepared this solution remains stable when stored at +4°C.
- **Control sera (+) and (-):**
Both controls are ready to use. Do not dilute
- **Conjugate:**
This reagent is ready to use. Do not dilute.
- **Serum diluent:**
This reagent is ready to use. Do not dilute.

VII. TEST PROCEDURE

1. Prior to starting the test, bring all reagents to room temperature (22-25°C).
2. Add 100 μ l of positive control to two wells of the plate, 100 μ l of the negative control to another 2 wells and 100 μ l of each of the dilutions of sera to be tested (prepared according to previous instructions) in the remaining wells of the plate. Seal the plate and incubate for **30 minutes at room temperature** (18-25°C).
3. Wash 3 times following the procedure previously described.
4. Add 100 μ l of the ready to use conjugate to each well. Seal the plate and incubate for **30 minutes at room temperature** (18-25° C).
5. Wash 5 times following the described procedure.
6. Add 100 μ l of the substrate solution to each well. Keep the plate at **room temperature for 15 minutes**. In order to speed up this process, it is advisable to use a multichannel pipette. Add 100 μ l of the stop solution to each well. We recommend adding this reagent following the same order in which the substrate was added
7. Read the OD of each well at 450 nm.

VIII. READING AND INTERPRETATION OF THE RESULTS:

Validation of the test:

The test is considered valid when the OD of the positive control at 450 nm minus the OD of the negative control is higher than 0.35 and the negative control OD at 450 nm is lower than 0.35.

Samples with a S/P higher or equal than 0.4 should be considered **positive** to PRRS antibodies.

Samples with a S/P lower than 0.4 should be considered **negative** to PRRS antibodies.

Interpretation of the results:

The S/P must be calculated as follows:

$$S/P = \frac{(\text{Sample OD}) - (\text{Negative C. OD})}{(\text{Positive C. OD}) - (\text{Negative C. OD})}$$

Negative S/P values must be considered with value 0.

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INMUNOLOGIA Y GENETICA APLICADA, S.A.