INgezim® RHDV1/2 DIF CROM

R.17.RHD.K.42

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INgezim® RHDV1/2 DIF CROM

R.17.RHD.K.42

INgezim® RHDV1/2 DIF CROM is a duplex immunochromatographic assay, which uses two antibodies specific to Rabbit Haemorrhagic Disease Virus (RHDV).

KIT FEATURES


APPLICATION

Detection and specific differentiation of Rabbit Haemorrhagic Disease Virus (RHDV and RHDV2) in rabbit liver homogenates and exudates.

TECHNICAL BASE

The diagnostic device consists of a conjugate pad that contains an antibody that recognizes the two RHDV variants, covalently bound to red latex particles; and a control protein bound to blue latex particles, which acts as an indicator of the correct development of the immunochromatography.


The membrane contains two test lines (T) and a control line (C): T1, where an antibody that recognizes the two RHDV variants is placed; T2, where a RHDV type 2 specific antibody is placed; and a control line, where an antibody specific to the control protein is placed. If the sample contains virus, it will react with the red latex particles bound to the antibody common to both variants, then, the formed complex will migrate through the membrane and will bind to the T1 line if the sample is RHDV type 1 positive or to both test lines (T1 and T2) if the sample is RHDV type 2 positive. In some cases, when the RHDV2 viral load is low, only the lower line (T2) may appear. In this way, depending on the number and position of red test lines that appear, the variant present in the sample can be determined. If the sample does not contain RHD virus, no colour will appear in any T line.

ASSAY VALIDATION

ANALYTICAL SENSITIVITY

In order to determine the analytical sensitivity of the test, serial dilutions of Virus Like Particles (VLPs), RHDV-VLPs and RHDV2-VLPs, were performed in running buffer and consequently analysed. As shown in the following figure, the limit of detection (LoD) of the assay was 0.08 µg/mL for RHDV-VLPs and 0.02 µg/mL for RHDV2-VLPs.

DIAGNOSTIC SENSITIVITY AND SPECIFICITY

A total of 96 liver samples were analysed, 23 corresponded to RHDV positive samples, 42 to RHDV2 positive samples and 31 to negative samples. Furthermore, a total of 40 liver exudates were analysed, 24 corresponding to RHDV2 positive samples and 16 to negative samples. Samples were previously classified as negative or positive to RHDV or RHDV2 by RT-PCR. Furthermore, liver samples were also analysed by a commercial ELISA for a common detection of RHDV and RHDV2 antigen (INgezim® RHDV DAS, Eurofins-Ingenasa, S.A).

Considering all the liver samples, results indicated >99.9% specificity and 93.8% sensitivity compared to RT-PCR. These results were compared with the results obtained when analysing the samples by the commercial ELISA described above, obtaining an excellent agreement (weighted kappa of 0.977%).

Regarding liver exudates, sensitivity reached 95.8% among RHDV2 samples, maintaining the specificity (>99,9%).

Detection	Liver samples		Liver exudates samples	
	Sensitivity % (95% CI)	Specificity % (95% CI)	Sensitivity % (95% CI)	Specificity % (95% CI)
RHDV	91.3 (71.9 – 98.7)	>99.9 (88.7-100.0)	ND	>99.9 (79.2-100.0)
RHDV2	95.5 (84.5 – 99.3)	>99.9 (88.7-100.0)	95.8 (78.8 – 99.3)	>99.9 (79.2-100.0)

COMPOSITION OF THE KIT

- Chromatographic devices
- Vials with diluent

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EUROFINS-INGENASA, S.A.