# INgezim® Maedi Screening

R.13.SLS.K.1

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## INgezim® Maedi Screening

R.13.SLS.K.1

**Ingezim® Maedi Screening** is an immunoenzymatic assay based on indirect ELISA technique, which uses a monoclonal antibody (MAB) specific to small ruminant immunoglobulins and 3 peptides specific to Maedi Visna virus (MVV) and to Caprine Arthritis-Encephalitis virus (CAEV).


### KIT FEATURES

#### APPLICATION

Detection of specific antibodies to Lentivirus in ovine and caprine serum samples.

#### TECHNICAL BASE

1. Plates supplied are coated with recombinant 3 peptides specific of Lentivirus (MVV and CAEV). On these wells, samples are added and incubated.
2. If samples contain specific antibodies to Lentivirus, they will bind to the antigen.
3. When the conjugate (monoclonal antibody specific to small ruminant IgG, conjugated with peroxidase, AcM-PO) is added, it will bind to the antibodies bound to the antigen. This binding is revealed by colorimetric reaction after addition of the substrate.

#### RESULTS INTERPRETATION

The assay two cut offs are used for the results interpretation. Samples are considered **Positive** if their IP value is higher than the positive cut off; **Negative** if their IP value is lower than the negative cut off and **Doubtful** if IP value is within the range of both cut offs.

### ASSAY VALIDATION

#### Diagnostic sensitivity and specificity

A panel of 84 sera from animals infected with genotypes A and B were analyzed. The results obtained indicated **95% sensitivity**.


On the other hand, a panel of 70 sera (genotype E infected and negative sera) were analyzed. The results obtained indicated a **specificity higher than 99.9%**.

#### Analytical specificity

It has been described in several references that some kind of cross-reaction may appear in assays for the detection of specific antibodies to M.visna in BTV vaccinated animals. In order to check the performance of the INgezim® Maedi Screening kit regarding this effect, different extractions at different days post vaccination of 20 animals vaccinated to BTV were analyzed. The results obtained indicated that INgezim® Maedi Screening does **not show cross-reaction** to this kind of antibodies.

#### Correspondence with other assays available in the market

44 ovine and caprine field sera were analysed by INgezim® Maedi Screening and two other commercial assays (1 and 2). The results obtained indicated **73% and 75% correspondence** respectively.

### KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (ABTS)
- Bottles with Stop Solution

Spanish registration n° 3295RD

EXPIRATION: 24 MONTHS. Stored at 2°C-8°C

**EUROFINS INGENASA, S.A.**

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28037 MADRID (SPAIN)

Phone: (+34)91 3680501

[www.ingenasa.com](http://www.ingenasa.com)

IT-73840  
IT-73780

9191.INGE 9175.ING2

Ed.280322