INgezim® West Nile IgM

R.14.WNV.K.2

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INgezim® West Nile IgM

R.14.WNV.K.2

INgezim® West Nile IgM is an immunoenzymatic assay based on the capture ELISA technique which uses a monoclonal antibody specific to equine IgM and a monoclonal antibody specific to West Nile Virus (WNV) E protein (Domain III).


KIT FEATURES

APPLICATION

Detection of West Nile IgM specific antibody in equine serum samples.

TECHNICAL BASE

1. Plates are supplied coated with a MAb specific of equine IgMs. Samples are added by duplicate to the wells and incubated.
2. If the samples contain IgMs (WNV specific or not), they will bind to the MAb coated.
3. Once the WNV antigen is added, the E protein will be captured by the specific IgM from the sample. A negative antigen is also used to determine unspecific reactions.
4. The WNV specific MAb-HRPO added later will bind to the antigen in case it has been captured by the specific IgM of the sample (infected or vaccinated animals). This binding is detected by adding a specific substrate.

RESULTS INTERPRETATION

The assay establishes two cut offs, which will classify the samples as **Positive** or **Negative**, depending on the value of the optical density of the sample in the assay, considering a range of optical densities close to the cut offs as **Doubtful** results.

ASSAY VALIDATION

VALIDATION WITH APHIS REFERENCE SERA

Two reference Sera from APHIS (USDA) were analysed:

- Equine serum infected with strain N.Y.99, positive by seroneutralization.
- Equine serum infected with strain N.Y.99, IgM positive.

Both sera were detected as positive.

ANALYTICAL SENSITIVITY

Different extractions made to vaccinated animal on days 0, 14, 21, 28, 42, 56, 70 and 85 post vaccination were analysed.

The assay was able to detect specific IgM to WNV from **day 14 post vaccination.**


RING TRIAL (Proficiency test, ANSES)

The INgezim® West Nile IgM assay participated in two Ring Tests organized by Anses (2013 and 2020). The results that were obtained are showed below:

- 2013: The conclusions of the final report regarding sensitivity were that INgezim® West Nile IgM kit is **more sensitive** than all the commercial assays available on the market that also participated in the Ring Test.
- 2020: According to the final report, INgezim® West Nile IgM is **more sensitive** than another commercial assay available on the market with which it was compared.

DIAGNOSTIC SENSITIVITY AND SPECIFICITY (ROC ANALYSIS)

Study 1: 14 equine sera infected by WNV, 1 SN positive reference serum, 1 WNV specific IgM positive serum, 3 equine WNV vaccinated sera positive by other commercial assays, 404 equine sera from WNV free areas and 28 field WNV infected sera from 3 outbreaks in Spain (2010-2012) and positive by another commercial ELISA were analysed. The sensitivity and specificity obtained were higher than 99.5%.

Study 2: 17 equine sera from a Linage 2 outbreak occurred in Italy in 2013 and previously classified as positive were analysed. INgezim® West Nile IgM detected 16 of them as positive.

KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Positive Antigen
- Vials with Negative Antigen
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with Stop Solution

Spanish registration nº 2863 RD
EXPIRATION: 18 MONTHS. Stored at 2°C-8°C

EUROFINS INGENASA, S.A.
C. de los Hermanos García Noblejas 39, 8º
28037 MADRID (SPAIN)
Phone: (+34)91 3680501
www.ingenasa.com