INgezim® West Nile Compac

R.10.WNV.K.3

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INgezim® West Nile Compac

R.10.WNV.K.3

INgezim® West Nile Compac is an immunoenzymatic assay based on a blocking ELISA technique, which uses a monoclonal antibody (MAb) specific to West Nile Virus (WNV) E protein.


KIT FEATURES

APPLICATION

Detection of specific antibodies to West Nile Virus (WNV) E protein in serum samples from horses and birds.

TECHNICAL BASE

- Plates are supplied coated with a semipurified extract of WNV proteins. Samples are added to the wells and incubated.
- If the samples contain specific antibodies to WNV, they will bind to the antigen.
- When the conjugate (monoclonal antibody specific to WNV E protein, conjugated with peroxidase, AcM-PO) is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the protein. In case the sample contains antibodies blocking the antigen (positive animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

RESULTS INTERPRETATION

The assay establishes two cut offs, which will classify the samples as **Positive** or **Negative**, depending on the IP value of the sample in the assay, considering a range of IP close to the cut offs as **Doubtful** results.

ASSAY VALIDATION

Evaluation with reference sera

The internal validation was performed with 3 OIE reference sera from laboratories of the national veterinary services USDA, APHIS. All of them showed the expected results.

Analytical sensitivity

Samples from an experimental infection of 6 partridges (isolate Spain/2007), 3 partridges (isolate Morocco /2003), 2 partridges as control and 2 rabbits (WNV NY99- crow) were analyzed. Extractions were made at days -1, 3, 6, 6, 10, 12, 14 and 20 or 0, 15, 24, 39 and 49 post-infection on partridges and rabbits, respectively.

The results obtained indicate that the assay detects specific antibodies of WNV since day 6 post infection.

Analytical Specificity

The assay does not show cross-reaction with specific antibodies to TBEV and SLEV.

Cross-reaction with specific antibodies to Japanese Encephalitis virus has been detected.

Regarding USUTU, cross-reaction has been detected in 1 of 13 samples analyzed.

RING TRIAL (Proficiency test ANSES 2013,2020)

According to the Final Report, INgezim® West Nile Compac showed **higher sensitivity and specificity** than another commercial assay available on the market with which it was compared.

Diagnostical sensitivity

Sera from 5 horses and 25 birds naturally infected with WNV (positive by seroneutralization (SN)) were analyzed.

Results obtained indicated that all of them were positive. The diagnostic sensitivity in this study was higher than **99.9%**.

Diagnostic specificity

86 sera from horses from Doñana, 247 sera from jackdaws, 254 sera from flamingos (all negative for SN) and 512 sera from horses collected during 2001, 2002 and 2004 in areas of Spain with no previous report of WNV seropositivity, were analyzed.

The diagnostic specificity was 96% for horses and 98% for birds. Hence, in general, the assay has **98.4% specificity**.

Comparative with other commercial kits

143 bird sera (20 positives and 123 negatives by SN) were analyzed and the results were compared with those obtained by another commercial assay available on the market. **The correspondence percentage was 91.6%.**

Compared to SN, **higher than 99.9% sensitivity and 86.2% specificity were observed.**

KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

Spanish registration nº 1815RD

EXPIRATION: 24 MONTHS. Stored at 2°C-8°C

EUROFINS INGENASA, S.A

C. de los Hermanos García Noblejas 39, 8º
28037 MADRID (SPAIN)
Phone: (+34)91 3680501
www.ingenasa.com