INgezim® ASFV CROM Ag

R.11.ASF.K.42

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INgezim® ASFV CROM Ag

R.11.ASF.K.42

INgezim® ASFV CROM Ag is a double antibody lateral flow assay based on the use of a monoclonal antibody (MAB) specific to VP72 of African Swine Fever Virus (ASFV).

KIT FEATURES


APPLICATION

Detection of African Swine Fever Virus (antigen) in blood samples of domestic pigs and wild boar.

TECHNICAL BASE

This diagnostic device consists of a plastic cassette containing a membrane and conjugated latex beads. On the membrane, two lines are printed: the test line (T) has a specific MAb to VP72 of ASFV and the control line (C) has a specific MAb to the control protein.

If a sample is positive, the virus binds to the black beads conjugated to the ASFV-specific MAb. The immune complex then migrates through the membrane by capillarity, and it is captured again by the ASFV-specific MAb absorbed on the test line, resulting in the appearance of a black line. The presence of the control line indicates the test is valid, the lateral flow assay has been performed correctly.


ASSAY VALIDATION

DIAGNOSTIC SENSITIVITY AND SPECIFICITY

To evaluate assay specificity, a total of 331 negative EDTA/heparin blood samples collected from domestic pigs and wild boar were evaluated. A 99.7 % specificity was obtained for the assay. In parallel, 14 EDTA blood samples characterised as negative by the reference technique (UPL-PCR) were analysed, obtaining a 100 % concordance.

To determine assay sensitivity, an external evaluation was carried out in collaboration with CISA-INIA, CSIC. A total of 67 EDTA blood experimental samples collected from domestic pigs and wild boar, characterised as positive by the reference technique (UPL-PCR), were evaluated. Positive samples were classified according to disease syndrome in acute, subacute, and chronic. An agreement with UPL of 61 % was found in the acute form of the disease, a 4 % in samples collected in the subacute clinical form, and a 10 % in the chronic clinical form.

ANALYTICAL SENSITIVITY

To determine the analytical sensitivity, blood negative samples spiked with recombinant VP72 or with inactivated virus (strain BA71) were evaluated. Serial dilutions in negative blood were made and the limit of detection was calculated.

INgezim® ASF CROM Ag was able to detect 38 ng/test of purified VP72 and down to 1/160 dilution of viral isolate.

COMPOSITION OF THE KIT

- Diagnostic devices
- Droppers with running buffer
- Graduated micropipettes

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EUROFINS-INGENASA, S.A.