INgezim[®] RHDV1/2 DIF CROM

Prod Ref: 17.RHD.K42

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I. TECHNICAL BASIS

The INgezim® RHDV1/2 DIF CROM test is a duplex lateral flow assay based on the use of coloured latex microspheres bound to the proteins of interest: red particles are covalently bound to an antibody that recognizes both RHDV variants and blue particles are coated with a control protein that indicates that the assay develops correctly.

The membrane contains two test lines (T) and one control line (C): T1, in which an antibody that recognizes both RHDV variants is immobilised; T2, that contains a RHDV2-specific antibody; and a control line (C), containing a control protein-specific antibody. If the virus is present in the sample, it will react with the red latex particles coated with the antibody for both variants and, subsequently, the complex will migrate across the membrane and bind to the protein contained in line T1 if RHDV type 1 or to both lines (T1 and T2) if it is RHDV type 2. In some cases, when the viral load of RHDV2 is low, only the bottom line (T2) can appear. In that way, according to the number and position of the red test lines that show up, the variant present in the sample can be determined. If the sample has no RHDV, no colour will appear in either T line.

II. STORAGE

All the reagents must be kept between +4°C and +25 °C.

III. NECESSARY COMPONENTS FOR PROCEDURE A NOT INCLUDED

- Weighing scale
- Centrifuge
- Petri dishes for sample processing
- Vortex
- Surgical scissors and / or 5 mL syringes

IV. PRECAUTIONS

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the kit.
5. Do not use components after expiration date and do not mix components from different batches.
6. Do not eat, drink, or smoke where specimens or kit reagents are being handle.

V. SAMPLES

The test has been designed for its use with samples of rabbit liver and liver exudates. Samples can be fresh or frozen at -20°C o -80°C. If using a frozen liver, it is important that it thaws completely before its extraction.

VI. TEST PROCEDURE

■ **Allow all kit components and samples to reach room temperature, before use.**

- If a liver sample is employed, the extraction can be carried out in a laboratory (part A of this section), or in field without any additional component (part B).
- If a liver exudate is employed, no sample extraction is required, go straight to part C of this section.

A. SAMPLE HOMOGENIZATION AND EXTRACTION IN LABORATORY

- Weigh one amount of liver (i.e. 1 g), place it in a Petri dish and cut it and squash it as much as possible using scissors and / or the syringe piston.
- Add the same volume (w/v) (i.e. 1 mL) of extraction buffer and continue squashing to form an homogenate.
- Using a pipette, take 0.5 mL of homogenate and place it in the supplied tube, in which 0.5 mL of extraction buffer have been dispensed previously.
- Vigorously stir for 30 seconds using a vortex and centrifuge 5 min at 4500 rpm.

- Take the supernatant and keep it at 4 °C if it is not going to be used immediately. For a longer period (more than 24 h), keep at -20°C.

B. SAMPLE EXTRACTION IN FIELD, WITHOUT ANY OTHER COMPONENT

- Take a small piece of liver and place it in the supplied tube (it should not surpass the 0.5 mL mark). If some scissors are available, cut the sample.
- Add the same volume of extraction buffer using a yellow pipette.
- Vigorously stir by hand for 30 seconds and leave for 5 min to allow phase separation.
- Take the supernatant and keep it at 4 °C if it is not going to be used immediately. For a longer period (more than 24 h), keep at -20°C.

C. ASSAY PROCEDURE

- Remove the diagnostic device from the bag and place it on a flat surface.
Add, using a transparent pipette, 20 µL of extracted supernatant or liver exudate onto the sample round window of the device. **Important:** If not enough supernatant is obtained for taking 20 µL, add one or two more drops of extraction buffer above it. After that, take the 20 µL needed.
- Add 3-4 drops of the supplied running buffer (wait between drops for the liquid to be absorbed).
- Read the results 10 minutes after the addition of the running buffer.

THE VOLUMES AND TIMES INDICATED MUST BE RESPECTED FOR THE TEST TO WORK CORRECTLY.

If any sample does not migrate correctly, it must be centrifuged again. If following the protocol B, it is recommended to leave the sample to rest for 1 hour, to facilitate phase separation.

1. Sample preparation:

- Liver
 - Homogenization and extraction in laboratory (A)
 - Extraction in field (B)
- Liver exudate

2. Add 20 µL of the extraction /exudate to the device.

3. Add 3 drops of running buffer and read the results after 10 min.

VII. INTERPRETATION OF RESULTS


Interpretation of results should be done as follows, 10 minutes after the addition of the sample:

POSITIVE result to RHDV1: (presence of RHDV1 in the sample)
A BLUE line appears in "C" zone (control) and a RED line appears in "T1" zone (test), the line situated in the middle.

POSITIVE result to RHDV2: (presence of RHDV2 in the sample)
A BLUE line appears in "C" zone (control) and 1) two RED lines appear in "T1" and "T2" zones (test) when viral load is high; 2) one RED line appears in "T2" zone (test) when the viral load is low.

NEGATIVE result: (absence of RHDV 1 or 2 in the sample). Only one BLUE line appears in the "C" zone (control).

INVALID result: Test is not valid if a BLUE line does not appear in position "C", either a "T" line appears or not in the test zone. Not valid results may be caused by damaged reagents or by an incorrect manipulation. In this case, the test should be repeated using a new strip.

It is advisable not to classify any sample as completely negative until the indicated time has elapsed since the analysis.

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INMUNOLOGIA Y GENETICA APLICADA, S.A.