# INgezim® FVR Compac

R.13.FVR.K3

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# INgezim® FVR Compac

R.13.FVR.K3

INgezim® FVR Compac is an immunoenzymatic assay based on a blocking ELISA technique, which uses a monoclonal antibody (MAb) specific to Rift Valley Fever Virus (RVFV) N protein and the recombinant N protein of the RVFV.


## KIT FEATURES

### APPLICATION

Detection of specific antibodies specific of Rift Valley Fever Virus N protein in multispecies assay for ruminants serum samples.

### TECHNICAL BASE

- Plates are supplied coated with antigen (recombinant N protein of RVFV). Samples are added to the wells and incubated.
- If the samples contain specific antibodies to RVFV, they will bind to the antigen.
- When the conjugate (monoclonal antibody specific to RVFV N protein, conjugated with peroxidase, AcM-PO) is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the antigen. In case the sample contains antibodies blocking the antigen (infected animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

### RESULTS INTERPRETATION

The assay establishes two cut off, which will classify the samples as **Positive** or **Negative**, depending on the Inhibition Percentage (IP) value of the sample in the assay, considering a range of inhibition percentages close to the cut off as **Doubtful** results.

## ASSAY VALIDATION

### DIAGNOSTIC SENSITIVITY


The following were analysed:

- 23 sheep experimentally infected with RVFV
- 8 previously sheep infected with RVFV.

The results obtained indicated 97% sensitivity.

### DIAGNOSTIC SPECIFICITY

1120 animals (sheep, goats, and cattle) from RVFV free areas were analysed. The results obtained indicated specificity values higher than 99%.

### KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (ABTS)
- Bottles with Stop Solution

Registration number in Spain: 3300 RD  
EXPIRATION: 15 MONTHS.  
Stored at 2°C-8°C

GOLD STANDARD DIAGNOSTICS MADRID, S.A.  
C/Hermanos García Noblejas 39, 8º  
28037 MADRID (SPAIN)  
Tel: (+34)91 3680501  
<http://www.goldstandarddiagnostics.com/>