

## INgezim MAEDI CONFIRMATION

R.13.SLC.K1

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

## INgezim MAEDI CONFIRMATON

R.13.SLC.K1

**INgezim MAEDI SCREENING** is an immunoenzymatic assay based in indirect ELISA technique, which uses a monoclonal antibody (MAb) specific of small ruminant immunoglobulins and two peptides: specific of Maedi Visna virus (MVV) and specific of Caprine arthritis-encephalitis virus (CAEV).

### TECHNICAL BASIS OF THE KIT

1. Two types of plates are supplied: coated with positive antigen which consists of 3 peptides specific of lentiviruses (MVV y CAEV) and coated with negative antigen.
2. Sera are added to both kind of plates and incubated for 1 hour at 37°C. If sera contain antibodies specific of MVV, they will bind to the positive antigen. In case of they contain unspecific antibodies, they will bind to both kind of antigens at the different plates.
3. After washing, a monoclonal antibody specific of ruminant IgG and conjugated with HRPO is added. The conjugate will recognize the ruminant IgG.
4. Finally, a substrate specific of HRPO is added to develop color in positive Wells after 10 min. of incubation. The subtraction of the values obtained in the positive and negative plates will determine if the simple is positive to Lentivirus or not.

### APLICACION

Confirmation of presence of antibodies specific of Small Ruminant lentivirus (MVV y CAEV) in ovine and caprine sera.

### INTERPRETATION OF THE RESULTS

The assay uses one cut off for results interpretation. Samples are considered **Positive** if their IP value is higher than cut off. Samples are considered **Negative** if their IP value is lower than the cut off.

### VALIDACIÓN

#### 3. Diagnostic sensitivity

A panel of 28 sera from infected herds and not controlled regarding genotype circulating were analyzed. These sera showed strong positive, weak positive and doubtful results by INgezim Maedi screening. In addition, 83 sera from animals infected with serotypes A (n=42) or B (n=41) were analyzed. The assay confirmed 86% of non-controlled animals and 100% of infected ones.

#### 4. Diagnostic specificity

A panel of 63 sera negative by other assays, either from free areas or from animals infected with serotype E were analyzed. All these sera result to be negative by INgezim Maedi Confirmation.

### COMPOSITION OF THE KIT

- Microtitration plates of 96 wells coated with positive antigen
- Microtitration plates of 96 wells coated with negative antigen
- Vial with Positive Control.
- Vial with Negative Control.
- Vial with Conjugate.
- Bottles with Washing Solution concentrated.
- Bottles with diluent
- Bottle with substrate (ABTS).
- Bottle with stop solution

PRODUCT MANUFACTURED BY INGENASA  
REGISTRATION NUMBER 3299 RD

Shelf life: **18 months**  
Store at 2°C-8°C

Ed.020217