INgezim PEA DAS

R.14.PEA.K2

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

INgezim PEA DAS

R.14.PEA.K2

INgezim PEA DAS is based on a Double Antibody Sandwich (DAS) ELISA technique, which uses specific antibodies to African Horse Sickness Virus (AHSV) VP7 protein.

TECHNICAL BASIS OF THE KIT

1. Plates are coated with specific antibodies for VP7 protein of AHSV. Samples are added and incubated.
2. If the samples contain the antigen, they will bind to the specific antibodies for the VP7 protein coating the plates.
3. When antibodies (labelled with biotin) specific to VP7 protein of AHSV are added, they will bind to the antigen previously bound to the antibodies coating the plate.
4. In a second amplification step, Streptavidin-Peroxidase is added, which will bind to the biotin. Binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APPLICATION

Detection of African Horse Sickness Virus (AHSV) in spleen samples



INTERPRETATION OF THE RESULTS

Two cut offs are used for the results interpretation. Samples are considered **Positive** if their OD value is higher than the positive cut off; **Negative** if their OD value is lower than the negative cut off and **Doubtful**, if show an OD value within both cut offs.

INTERNAL VALIDATION


1. Use of culture samples of AHSV

The sensitivity of the assay is **15ng or 5×10^4 UFP/ ml** for purified AHSV and supernatant of AHSV infected cells samples.

2. Field trials

A panel of 91 equine spleen samples obtained from areas in which an outbreak appeared and from AHSV free areas were analyzed. Samples were previously classified by virus isolation. These results indicated that the correlation between both techniques was **93%**.

EXTERNAL VALIDATION

The assay has been evaluated in the Animal Health Laboratory of Maisons-Alfort concerning analytical sensitivity. The conclusions of the study indicated that INGEZIM PEA DAS is suitable for the identification of the nine serotypes of AHSV, provided that the viral load in tissues or culture supernatants exceed 10^3 TCID₅₀/ml.

COMPOSITION OF THE KIT

- Microtitration plate of 96 wells
- Vial with Positive Control
- Vial with Negative Control
- Vial with Biotin Conjugate
- Vial with Streptavidine-HRPO
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with substrate (ABTS)
- Bottle with stop solution
- Bottle with Substrate Buffer

PRODUCT MANUFACTURED BY INGENASA
REGISTERED IN SPAIN 0427-RD

SHELF LIFE: **18 months**
Stored at 2°C-8°C

Ed.250117