Avian metapneumovirus (APV)

## Anticuerpo Monoclonal específico de P de APV P APV specific Monoclonal Antibody

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*Avian pneumovirus (APV) causes a respiratory disease in turkeys (turkey rhinotracheitis) and others avian species (Swollen head syndrome, in chickens). It belongs to the genus Metapneumovirus. In turkeys, the disease is characterized by rales and nasal discharge, conjunctivitis, swelling of the infraorbital and submandibular sinuses. Mortality isn't higher than 2-5% although morbidity can reach up to 10%. It can decrease the eggs production in 10-30% with the consequent economic loss. According nucleotides sequence of the union G glycoprotein, 2 groups were originally described, A and B. The first American isolates were genetically and antigenically different so a new group was described, C. Finally, French isolates obtained in the 80's were different to the other 3 and a fourth group, D, was defined. APV is a simple RNA virus, with negative polarity and an envelope, which codifies for 8 structural proteins. The genome is 13.3 Kb. Glycoproteins G and F are the main neutralizing antigens and the N protein is the most conserved between different subtypes.*

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*The hybridome which produces the monoclonal antibody has been obtained by the fusion of lymphocytes from Balb/c mice's spleen with myeloma X63/Ag8653 cells. The IgG has been purified by HPLC presenting a purity of 99%, its isotype is IgG<sub>1</sub> and it is specific to APV's P protein.*

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
*Detection of APV by Immunohistochemistry (IHC) technique.*

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
*Detection of APV by indirect immunofluorescence technique (IFI).*

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*Nasal turbinates of the APV-infected chickens are fixed in neutral buffered formalin for 24-92 hr, and then embedded in paraffin by standard procedures. The slides obtained are deparaffined, and rehydrated through graded alcohols. The endogenous peroxidase activity is blocked by incubation with hydrogen peroxide and the slides are then immersed in 10 mM citrate buffer (pH 6) in bathwater at 98° C for 20 min. Incubate with 20% negative goat serum solution in 0.1 M Tris- buffered saline for 1 hour at room temperature, followed by another incubation with the mouse monoclonal antibody at 4° C overnight (See recommended use conditions for dilution. Note that the supernatant is supplied ready to use). As a secondary antibody, biotinylated goat anti-mouse at 1:200 is used, followed by an ABC complex diluted 1:100 in TBS, applied for 1 hour at room temperature. Sections are finally incubated in TBS with diaminobenzidine (DAB)-hydrogen peroxide solution for 10 min., counterstain with Harris's haematoxylin, dehydrate and cover with a cover slip.*

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