
INgezim CORONA CANINO

R.15.CCV.K1

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

INgezim CORONA CANINO

R.15.CCV.K1

INgezim CORONA CANINO is based on an indirect enzymatic immunoassay (indirect ELISA), which uses a specific monoclonal antibody (MAb), specific to canine immunoglobulins (Igs) and a semipurified extract of Canine Coronavirus (CCV) as an antigen.

TECHNICAL BASIS OF THE KIT

1. Plates are coated with CCV antigen. Serum samples are added and incubated.
2. If the samples contain specific antibodies to CCV, they will bind to the antigen.
3. When a MAb-PO specific to canine immunoglobulins is added, it will bind to the Igs previously bound to the antigen. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APPLICATION

Detection of specific antibodies to CCV in canine sera samples.

INTERPRETATION OF THE RESULTS

One cut off is used for the results interpretation: OD values higher than the cut off are considered **Positive**; OD values lower than the cut off are considered **Negative**.


VALIDATION

1. Sensitivity

In this study 16 sera from CCV vaccinated animals and 12 from CCV infected animals were analyzed being all of them detected as positive by INgezim CORONA CANINO.

2. Specificity

In this study 14 sera from CCV negative animals were analyzed being all of them detected as negative by INgezim CORONA CANINO.

COMPOSITION OF THE KIT

- 96 wells microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with conjugate
- Bottle with washing solution
- Bottle with diluent
- Bottle with substrate
- Bottle with stop solution

PRODUCT MANUFACTURED BY INGENASA

Shelf life: **18 months**
Store at: 2°C-8°C

Ed. 020217