

## INgezim FeLV DAS

R.16.FLV.K2

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

## INgezim FeLV DAS

R.16.FLV.K2

**INgezim FeLV DAS** is an enzymatic assay based on a Double Antibody Sandwich ELISA technique, which uses monoclonal antibodies (MAb) specific of Feline Leukaemia Virus (FLV) p27 protein.

### TECHNICAL BASIS OF THE KIT

1. Plates are supplied coated with a monoclonal antibody specific of FLV protein p27. Samples are added and incubated.
2. If samples contain FLV, they will bind to the specific monoclonal antibody.
3. After incubation, a MAb-HRPO specific of FLV p27 protein is added. The mix will bind to the protein if there are antigen in the sample
4. This union is detected by a colorimetric reaction after the addition of the substrate specific of HRPO.

### APLICACION

Detection of FLV antigen p27 in sera and plasma


### INTERPRETATION OF THE RESULTS

One cut off is used for results interpretation. Samples will be **Positive** if their OD value is higher than the cut off. Samples will be **Negative** if their OD value is lower than the cut off.

### VALIDATION

#### 1. Analytical sensitivity.

In order to determine the analytical sensitivity of INgezim FeLV DAS 3 samples were analyzed: viral antigen, recombinant protein (gag gen) and negative antigen. Results obtained indicated that the assay is able to detect 10ng of antigen in both kind of positive samples.

#### 2. Comparative with other commercial kits

55 feline sera previously classified by another commercial kit based on ELISA technique were analyzed. In this study, correspondence between both assays was 100%.

Ingezim FeLV DAS				
commercial kit		POS	NEG	TOTAL
	POS	11	0	11
	NEG	0	44	44
	TOTAL	11	44	55

### COMPOSITION OF THE KIT

- 96 wells microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with conjugate
- Bottle with washing solution
- Bottle with diluent
- Bottle with substrate
- Bottle with stop solution

PRODUCT MANUFACTURED BY INGENASA

Shelf life: **15 months**  
Store at: 2°C-8°C

Ed. 020217