

## INgezim RHDV

R.17.RHD.K1

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

## INgezim RABBIT

R.17.RHD.K1

**INgezim RABBIT** is based on an indirect ELISA technique, which uses Protein A conjugated with Peroxidase, and a recombinant antigen (VP1 protein of RHDV).

### TECHNICAL BASIS OF THE KIT

1. Plates are coated with RHDV antigen (recombinant protein VP1). Serum samples are added and incubated.
2. If the samples contain specific antibodies to RHDV, they will bind to the antigen.
3. When the Protein A-PO is added, it will bind to the Igs previously bound to the antigen. This binding is detected by the development of a colorimetric reaction after the addition of the substrate.

### APPLICATION

Detection and/or titration of specific antibodies to RHDV in rabbit sera samples


### INTERPRETATION OF THE RESULTS

There is one Cut Off for the results interpretation. Samples with an OD higher than Cut Off value are **Positive** and with OD lower than Cut Off value are **Negative**. The titre of the sample will be the last dilution of the sample showing an OD higher than the Cut Off value.

## VALIDATION


### ANALYTICAL SENSITIVITY

From studies made with animals intramuscularly vaccinated and challenged at day 36 post vaccination, it can be concluded that INGEZIM RHDV is able to detect specific antibodies to VP60 since day 14 p.v. depending on the vaccine and the dose used. (0.5-3-10-25-50µg; commercial vaccine).

### CORRELATION WITH THE REFERENCE TECHNIQUE IHA (INHIBITION OF HAEMAGLUTINATION).

- **Study 1:** Titres obtained by both techniques at different days p.v. from 4 orally vaccinated animals were compared. In this study it was observed that INGEZIM RHDV was able to detect specific antibodies in all of the studied animals since day 21 p.v. while IHA detected these antibodies in 3 animals.

- **Study 2:** A panel of 150 field sera was analyzed by INGEZIM RHDV and IHA. The results obtained indicated that INGEZIM RHDV showed 99% sensitivity in relation to IHA and 86% specificity, which was probably due to the fact that the ELISA is able to detect specific antibodies earlier than IHA does, as it has been demonstrated in other studies.

### COMPOSITION OF THE KIT

- Microtitration plates of 96 well coated with antigen
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate (-20°C)
- Bottle with Washing Solution
- Bottle with Diluent
- Bottle with Stop Solution
- Bottle with Substrate (ABTS)
- Bottle with Substrate Buffer

PRODUCT MANUFACTURED BY INGENASA  
SPANISH REGISTRATION NO: 3198 RD

SHELF LIFE: **18 months**  
Stored at 2°C-8°C

Ed.250117