

INgezim PPV R.11.PPV.K1

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

CADUCIDAD: **18 meses**
Conservado a 2°C-8°C

INgezim PPV R.11.PPV.K1

INgezim PPV is based on an indirect ELISA technique, which uses a monoclonal antibody (MAb) specific to porcine immunoglobulines (Igs), and recombinant antigen (VP2 protein of PPV).

TECHNICAL BASIS OF THE KIT

1. Plates are coated with PPV antigen (recombinant capsides of VP2). Serum samples are added and incubated.
2. If the samples contain antibodies specific to PPV, they will bind to the antigen.
3. When a MAb-PO specific of porcine immunoglobulines is added, it will bind to the Igs bound to the antigen. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APPLICATION

Detection and /or titration of specific antibodies to PPV in porcine sera samples.

INTERPRETATION OF THE RESULTS

There is one Cut Off for the results interpretation. Samples with an OD higher than the Cut Off value are positive and samples with an OD lower than the Cut Off value are negative. The titer of the sample will be the last dilution of the sample showing an OD higher than the Cut Off value

SENSITIVITY AND SPECIFICITY


1. Correspondence with the reference technique IHA (Inhibition of haemagglutination)

238 field sera were analysed, to determine the sensitivity and the specificity of INGEZIM PPV ELISA technique in relation to the IHA technique

SENSITIVITY	SPECIFICITY
99.5%	100%


3. Sensitivity using experimentally vaccinated sera.

Experimentally vaccinated animals were used and the antibodies detection profile was obtained. In the same study, a comparison with the IHA results was made. The results are represented in the following figure.

2. Analytical sensitivity

On the second week post infection, the assay can detect titres of antibodies around 1/200.

4. Specificity using SPF animals

In order to determine the specificity of the assay, a set of 188 SPF negative animals was used. The results obtained indicated 100% specificity.

5. One well titration of sera

It is possible to obtain an approximated titre of the sera based on its OD value, by using a formula which corresponds with an exponential curve.

COMPOSITION OF THE KIT

- Microtitration plates of 96 wells coated with Ag
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with Diluent
- Bottle with Substrate (ABTS)
- Bottle with Substrate Buffer
- Bottle with Stop Solution

PRODUCT MANUFACTURED BY INGENASA
REGISTRATION NUMBER: 1129 RD

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