

INgezim PPV DAS

R.11.PPV.K2

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

INgezim PPV DAS

R.11.PPV.K2

INgezim PPV DAS is an enzymatic assay based on a Double Antibody Sandwich ELISA technique, which uses monoclonal antibodies (MAb) specific to Porcine Parvovirus (PPV) VP2 protein.

TECHNICAL BASIS OF THE KIT

1. Plates are coated with a monoclonal antibody specific for VP2 protein of PPV. On these wells, samples are added and incubated.
2. If the samples contain the antigen, it will bind to the specific antibodies for the VP2 protein which are coating the plate.
3. When a MAb-BIOTIN specific for VP2 protein of PPV is added, it will bind to the antigen bound to the antibodies coating the plate.
4. In a second amplification step, Streptavidin-Peroxidase is added, which will bind to the biotin. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APPLICATION

Detection of Porcine Parvovirus (PPV) in foetal tissue samples.

INTERPRETATION OF THE RESULTS

One cut off is used for the results interpretation. Samples will be **Positive** if their OD value is higher than the cut off. Samples will be **Negative** if their OD value is lower than the cut off.

VALIDATION

SENSITIVITY AND SPECIFICITY

1. Correspondence with the reference technique HA (haemagglutination)


A set of 189 samples were analyzed. 120 samples belonged to foetus from vaccinated gilts which were infected with PPV. The other 69 samples belonged to foetus from non vaccinated and PPV infected gilts. The results indicated 97% and 98% specificity and sensitivity respectively in relation to HA.

2. Sensitivity of the assay

To determine the sensitivity of the assay, different dilutions of a standard of purified PPV VLPs were titrated. The concentration of the standard was determined by amino acid analysis. The results obtained indicated **that the assay is able to detect 1ng of standard of PPV VLPs.**

3. Samples quantification

The ELISA DAS can be used to determine precisely the amount of PPV VLPs present in a sample, by using a formula which corresponds with a linear regression curve.

COMPOSITION OF THE KIT

- Microtitration plate of 96 wells
- Vial with Positive Control
- Vial with Negative Control
- Vial with Biotin Conjugate
- Vial with Streptavidine-HRPO
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with substrate (ABTS)
- Bottle with Substrate Buffer
- Bottle with stop solution

PRODUCT MANUFACTURED BY INGENASA

SHELF LIFE: **18 months**
Stored at 2°C-8°C

Ed.020217