

## INgezim PPV Compac

R.11.PPV.K3

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

## INgezim PPV Compac

R.11.PPV.K3

**INgezim PPV COMPAC** is an enzymatic assay based on a blocking ELISA technique which uses a monoclonal antibody (MAb) specific to porcine parvovirus (PPV) VP2 protein, and a recombinant capsid of VP2.

### TECHNICAL BASIS OF THE KIT

1. Plates are coated with PPV antigen (recombinant capsides of VP2).
2. Serum samples are added and incubated.
3. If the samples contain antibodies specific to PPV only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the protein. In case the sample contains antibodies blocking the antigen (infected animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

### APPLICATION

Detection of specific antibodies to porcine Parvovirus, in porcine sera samples.

### INTERPRETATION OF RESULTS

Two Blocking Percentage Values are used for the results interpretation. Samples with a Blocking Percentage higher than 30% are considered **Positive** and samples with a Blocking Percentage lower than 25% are considered **Negative**. Samples with a Blocking Percentage between both values must be considered **Doubtful**.

## VALIDATION

### SENSITIVITY

1. Sensitivity with experimental sera  
A set of sera from 16 experimentally vaccinated animals and 10 not vaccinated animals was used. After vaccination, 12 different extractions of each animal were taken and the presence of antibodies was determined by INgezim PPV COMPAC and IHA. At day 84, all animals received a challenge and blood sampling was taken once a week during 6 weeks more. The results obtained indicated that INgezim PPV COMPAC was able to detect specific antibodies to PPV between days 7 and 15 after the antigen contact.

### SPECIFICITY

1. Specificity using negative controlled animals

A set of 153 sera from negative controlled animals was used. These animals had never been in contact with the virus (SPF). The results obtained indicated that the specificity of the assay was 100%.

### CORRESPONDENCE WITH THE REFERENCE TECHNIQUE IHA

A set of 317 sera from experimentally vaccinated and not vaccinated animals were analysed by INgezim PPV COMPAC and the results were compared with those obtained with the Inhibition of haemagglutination (IHA) assay. The results obtained indicated that the correspondence between INgezim PPV COMPAC and IHA was **99.7%**.

### COMPOSITION OF THE KIT

- Microtitration plates of 96 wells.
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with stop solution
- Bottle with substrate (TMB) ready to use

PRODUCT MANUFACTURED BY INGENASA

SHELF LIFE: **18 months**  
Stored at 2°C-8°C

Ed.020217