Anticuerpo Monoclonal específico de N de PRRSV N PRRSV specific Monoclonal Antibody

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The Porcine Reproductive and Respiratory Syndrome (PRRS) is a disease caused by an Arterivirus (PRRSV) widely spread in pigs all over the world. There are two clear different serotypes: the European and the American. The PRRSV is a single strain RNA virus with positive polarity and a size of 15kb. The genome has 8 Open Reading Frames (ORF) which codify for seven viral proteins. ORF 1_a and

1_b codify for the viral polymerase and the rest for the viral structural proteins. The virus has three main proteins: E or gp50, M and N which correspond respectively with the expression products of the ORF

5, 6 and 7. The N protein is the most immunogenic one but the specific antibodies against it can't be detected until day 14 post infection, and are not in any case neutralizing.

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
The hybridome which produces the monoclonal antibody has been obtained by the fusion of lymphocytes from Balb/c mice's spleen with myeloma X63/Ag8653 cells. The IgG has been purified by HPLC, showing a purity of 99%. Its isotype is IgG_{2a} and It is specific to the N protein to PRRSV.

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
Detection of PRRSV by Immunohistochemistry (IHC) technique.

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Samples are fixed in 10% saline formalin for 20-24 hours and processed to be included in paraffin. Obtain 3µm sections and transfer to slides, previously treated with Vectabond. Dry at 37°C over night. Samples are analyzed with or without treatment with protease. Antibody is used diluted in rabbit serum at 1% in PBS (See recommended use conditions for dilution. Note that the supernatant is supplied ready to use) Samples are incubated in wet chamber at 4°C for 18 hours. As a second antibody a biotinylated rabbit antimouse immunoglobulin is used. The reaction is amplified with streptavidine- peroxidase and diaminobenzidine is used as a chromogen.

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