
INGEZIM PRRS 2.0

R.11.PR2.K1

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
^aTatjana Sattler et al. BMC Veterinary Research 2014, 10:300.; ^bTatjana Sattler et al. Porcine Health Management (2015) 1:19

Estudio 2: 2014 International PTS for European and American PRRSV strains specific antibody detection (GD Animal Health, Deventer, The Netherlands). Participaron 87 laboratorios con ensayos diferentes y se analizaron 8 muestras. Los resultados se muestran en la figura:

2018 International PTS for European and American PRRSV strains specific antibody detection (GD Animal Health, Deventer, The Netherlands). 100% de correspondencia con los resultados esperados.

		IDEXX X3			Acuerdo (κ)	Coeficiente correlación
		POS	NEG	TOT		
INGEZIM	POS	112	8	120	0,89	0,90
	NEG	0	50	50		

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

INgezim PRRS 2.0

R.11.PR2.K1

INgezim PRRS 2.0 is based on an indirect ELISA technique, which uses a monoclonal antibody (MAb) specific for porcine immunoglobulins (IgG), and both recombinant N proteins of American and European strains of PRRS as antigens.

TECHNICAL BASIS OF THE KIT

1. Plates are supplied coated with PRRSV antigen (recombinant N protein from European and American strains). Serum samples are added and incubated.
2. If the samples contain specific antibodies to PRRSV, they will bind to the antigen.
3. When a MAb-PO specific of porcine immunoglobulins is added, it will bind to the IgG previously bound to the antigen. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APPLICATION

Detection of specific IgG to PRRSV in porcine sera samples. The assay has been designed for the detection of IgG specific for both, American and European strains of PRRSV.

INTERPRETATION OF THE RESULTS



The assay establishes one Cut Off value measured in S/P: Samples with an OD higher than the Cut Off value are considered positive and samples with an OD lower than the Cut Off value are considered negative.

ANALYTICAL SENSITIVITY

A. Vaccinated animals

- **European strain:** (Barcelona University). 192 animals were analyzed by INgezim PRRS 2.0. These animals were experimentally vaccinated with an European vaccine and bled at days 0, 7, 14, 21. **97% were detected positive at day 14 p.v.**
- **American strain:** (INgezim PRRS 2.0 (■) IDEXX X3 (■)) 3 experimentally vaccinated animals bled at different days p.v. (Ingelvac PRRSV) **Analytical sensitivity for American strain since day 12 p.i.**

- **B. Infected animals (European strain):** (INgezim PRRS 2.0 (■) IDEXX X3 (■)). 7 animals were bled at different days p.i. **Analytical sensitivity for European strain since day 11 p.i.**

DIAGNOSTIC SENSITIVITY AND SPECIFICITY

1615 American and European field sera were analyzed by INgezim PRRS 2.0 and by IDEXX X3. Sensitivity and specificity in comparison with IDEXX X3 were:

	EUROPEAN STRAIN	AMERICAN STRAIN
SENSIBILITY	99,3	91
SPECIFICITY	98,2	96,5


EXTERNAL VALIDATION

Study 1^a (detection of Ab specific of European strains): The following sera were analyzed by INgezim PRRS 2.0 and by IDEXX X3 as "gold standard": 3 pigs vaccinated with an attenuated vaccine genotype 2, 245 pigs from positive herds (group II), 309 pigs from negative herds (group III), 256 from expected positive herds without any history (group IV) and 92 wild boars (group V). Results obtained indicated:

- Animals vaccinated with genotype 2 vaccine: Coefficient κ respect to IDEXX X3 $>0,8$; detecting Ab in 1 pig at day 9 p.v. and in the rest at day 12 p.v.

		Group used
Sensibility	96%	I-II (n=266)
Specificity	99%	III (n=309)
Predictive positive value	86%	I-II-III-IV (n=828)
Predictive negative value	98,7%	I-II-III-IV (n=828)
Accuracy	95,3%	I-II-III-IV (n=828)
Coefficient Kappa (κ)	0,88	I-II-III-IV (n=828)

Study 3: 2014 International PTS 2014 for European and American PRRSV strains specific antibody detection (GD Animal Health, Deventer, The Netherlands). The study included 87 different laboratories using different ELISAs and analyzing 12 different samples. Results are shown in the picture:

2018 International PTS for European and American PRRSV strains specific antibody detection (GD Animal Health, Deventer, The Netherlands). 100% de correspondencia with the expected results.

	IDEXX X3			Agreement (κ)	Correlation Coefficient
	POS	NEG	TOT		
INGEZIM	POS	112	8	0,89	0,90
	NEG	0	50		

Study 3^b (detection of Ab specific of Asian strain): 10 piglets were injected with an attenuated PRRSV type 2 vaccine "Kyoto Biken" (PIGWIN PRRSV2) and other 10 were used as control. These 20 animals were challenged at the aged of 10 weeks with an HP PRRSV field strain (Vietnam_PRRSV_AGES/568-30FC/13). **Results obtained indicated a correlation coefficient of 0.90 with $\kappa=0.89$ of agreement to detect Ab specific of the Asian strain in comparison with IDEXX X3.**

^aTatjana Sattler et al. BMC Veterinary Research 2014, 10:300. ^bTatjana Sattler et al. Porcine Health Management (2015) 1:19

COMPOSITION OF THE KIT

- Microtitration plates of 96 wells
- Vials with Positive Control, ready to use
- Vials with Negative Control, ready to use
- Vials with Peroxidase conjugate ready to use
- Bottle with Washing Solution
- Bottle with Diluent for sera
- Bottle with stop solution
- Bottle with substrate (TMB)

MANUFACTURED BY INGENASA

SHELF LIFE: **15 months**
Stored at 2°C-8°C

Ed.080518