INGEZIM FVR Compac

Prod Ref: 13.FVR.K3

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I. TECHNICAL BASIS

This kit is based on a blocking enzymatic immunoassay (Blocking Elisa). The technique is briefly described below:

The antigen is fixed in a solid support (polystyrene plate). After incubating the sera samples, a RVFV specific monoclonal antibody (Mab peroxidase conjugated) is added. If the sample contains antibodies specific of the virus, they will not allow the binding of the labelled Mab to the antigen whereas if it does not contain specific antibodies, the Mab will bind to the antigen which is coating the plate.

After washing the plate to eliminate all non-fixed material, the presence or absence of labelled Mabs can

be detected by adding the substrate which, in presence of the peroxidase, will develop a colorimetric reaction.

If there has been colour development, this means that the conjugate has bound to the antigen, being therefore the sample negative. On the other hand, if the sample contains specific antibodies to RVFV, they will block the binding of the conjugate and there will no colour development will be detected.

The antigen used in this kit is the N recombinant protein from the RVFV. This confers many advantages to the assay. Being produced recombinantly, it therefore totally lacks of infectivity. Plus the N protein is one of the mayor proteins from RVFV, as well as the most antigenic.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix reagents or use instructions from different kits.
4. Avoid any contamination of the kit's reagents.
5. Do not use components after expiry dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking while handling reagents or samples.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Systematically include a positive control and a negative serum each time the assay is run.
10. Substrate must be handled with care, as it is very sensitive to light and contamination.
11. Stop solution can precipitate at refrigeration. Make sure there is no precipitation before its use.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored at +4°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps can be done using an automatic plate washer or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the instructions below:

- Throw away the content of the plate by briskly turning the plate over, in order to avoid the possible exchange of contents from one well to another
- Dispense a volume of 300 µl of washing solution in each well.
- Shake the plate gently, avoiding cross-contamination between wells.
- Briskly turn the plate over to empty the wells.
- Repeat the process as many times as indicated in the kit's instructions.
- Prior to emptying the content of the plate in the last washing step, verify that the next reagent to be added is ready to be used. Do not let the plate dry longer than strictly necessary to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use.
- After the last washing step tap the plate upside down on absorbent filter paper to remove any remaining washing solution

V. PREPARATION OF REAGENTS

○ **Washing solution:**

Dilute one part of the concentrated washing solution provided in the kit with 24 parts of distilled or deionised water (i.e. 40 ml of concentrate and 960 ml of water). Once prepared, this solution remains stable when stored at +4°C.

○ **Preparation of Controls (+) y (-):**

Both controls are ready to use and do not require any further preparation. Dispense 100 µl of each one.

○ **Preparation of conjugate:**

This reagent is ready to use. Do not dilute

VI. PREPARATION OF SAMPLES:

Sera samples must be assayed at 1/5 dilution in the diluent provided. This dilution can be done directly in the plate by adding 80 µl of the diluent and then adding 20 µl of sample to each well.

VII. TEST PROCEDURE

1. Prior to starting the test, bring all reagents to room temperature (22-25°C).
2. Add 80 µl of the diluent and then 20µl of each sample in every well to be assayed. Add 100µl of positive and negative controls to their respective wells. We recommend running controls in duplicate. **Seal the plate and incubate 45 minutes at 22-25°C.**
3. Wash 4 times following the procedure previously described.
4. Add 100 µl of the ready to use conjugate to each well. **Seal the plate and incubate for 30 minutes at 22-25°C.**
5. Wash 4 times following the procedure previously described.
6. Add 100 µl of the substrate solution to each well. Keep the plate **at room temperature for 10 min.**
7. Add 100 µl of the stop solution to each well.
8. Read the OD of each well at 405 nm.

VIII. READING AND RESULT INTERPRETATION

When running samples in duplicate, the OD value for each sample should be calculated as the arithmetic mean of both values. In the same way, the arithmetic mean of the values obtained in the two wells of the positive control and the two wells of the negative control must be calculated.

A. VALIDATION OF THE RESULTS

The test is considered valid when:

- The OD of the Negative Control (NC) is higher than 0.8.
- The OD of the Positive Control (PC) is lower than 0.25.

B. INTERPRETATION OF THE RESULTS

Calculate the inhibition percentage (IP) of each sample as follows:

$$IP \text{ of sample} = 100 - [(sample \text{ OD} / NC \text{ OD}) \times 100]$$

Samples will be considered **POSITIVE** (i.e. the sample contains specific antibodies to RVFV), when the IP is $\geq 45\%$.

Samples will be considered **NEGATIVE** (i.e. the sample contains specific antibodies to RVFV) when the IP is $\leq 40\%$

Samples with an IP between 40 and 45% should be considered as **DOUBTFUL and be retested by a seroneutralization (SN) assay.**

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by: