INGEZIM BTV DR

Prod Ref: 12.BTV.K0

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I. TECHNICAL BASIS

The kit has been designed to detect antibodies specific for BTV in sheep, goats and cattle being able to detect a very low titers of antibodies in sera of infected and vaccinated animals.

The INGEZIM BTV DR kit is based on a novel immunoenzymatic assay called double recognition ELISA which is described below:

Plates are coated with VP7 protein of BTV. After adding the sample to the well, if it contains BTV specific antibodies, they will bind to the antigen coating the plate while if the sample does not

contain specific antibodies they will not bind to the antigen.

When VP7 protein conjugated with peroxidase is added and sera samples contain VP7 specific antibodies, they will catch the labelled VP7. Presence or absence of labelled VP7 will be detected by addition of substrate (TMB) which, in presence of the peroxidase, will develop a colorimetric reaction.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
10. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
11. Substrate must be handled with care, it is very sensible to light and contamination.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionised water. Once prepared, this solution remains stable between +2°C and +8°C.

- **Positive and Negative controls:**

Controls are ready to use and must be tested adding 100µl / well.

- **Preparation of the conjugate:**

Conjugate is ready to use and must be used adding 100µl/well

VI. PREPARATION OF SAMPLES:

Sera samples must be assayed diluted 1/2 in diluent. This dilution could be made directly on the assay plate by adding 50 µl of diluent and 50 µl of sample to each well.

VII. TEST PROCEDURE

1. All reagents must be allowed to come to room temperature before use.
2. Addition of samples and controls:
 - **Samples:** 50 µl of diluent to each well. Add 50 µl of each sample to be assayed.
 - **Controls:** Add 100µl of positive control to two wells and 100 µl of negative control to other two wells. We recommend running controls in duplicate.

Seal the plate and incubate 1 hour at 37°C.
3. Wash 6 times following the described procedure.
4. Add 100 µl of conjugate ready to use to each well. **Seal the plate and incubate for 1 hour at 37°C.**
5. Wash 6 times following the described procedure.
6. Add 100 µl of substrate, to each well. Keep the plate for **15 min at room temperature.**
7. Add 100 µl of stop solution to each well.
8. Read the OD of each well at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

Determine the mean absorbance of positive and negative controls and samples, in case of assayed the last ones in duplicate.

1. VALIDATION OF THE RESULTS

Positive Control > 0.8
Negative Control < Cut off

2. CUT OFF CALCULATION

Cut Off = 0.15 x Positive control

3. RESULTS INTERPRETATION

Samples will be considered **POSITIVE** (there are antibodies specific of BTV), if the OD value at 450nm is higher than the cut off (15 % of positive control).

Samples will be considered **NEGATIVE** (there are no antibodies specific for BTV in the sample) if the OD value at 450nm is equal or lower than the positive cut off.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by: