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I. TECHNICAL BASIS

The kit has been designed for detection and quantification of BTV's VP7 protein.

The INGEZIM BTV DAS kit is based on the double antibody sandwich immunoassay technique it described below:

The solid phases are plates coated with a monoclonal antibody (MAb) specific of BTV's VP7 protein. After adding the sample to the well, if it contains virus, they will bind to the MAb absorbed to the plate .If the sample

does not contain virus they will not bind to the MAb. After washing to eliminate the not adsorbed material, a MAb specific to VP7 conjugated with peroxidase is added. If there was virus in the sample, the conjugate will bind to the antigen. This union is detected after adding a substrate able to develop colour in presence of peroxidase.

This ELISA can detect up to 3 ng/ml of VP7.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each sample.
9. For each utilisation of the Kit, control positive and negative sample must be tested in a systematic way.
10. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
11. Substrate must be handle with care, it is very sensible to light and contamination.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF REAGENTS

- **Washing solution:**
Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionised water (i.e. 40 ml of concentrate and 960 ml of water). Once prepared, this solution remains stable between +2°C and +8°C.
- **Diluent :**
Dilute 1 part of 5x concentrated diluent with 4 parts of distilled or deionised water.
- **Preparation of conjugate:**
It is ready to use. Do not dilute.
- **Preparation of standards and negative control.**
They are ready to use. Do not dilute.

VI. PREPARATION OF SAMPLES:

Samples must be assayed at two fold dilutions from 1/25 in the diluent supplied. It is recommended to make several dilutions because the absorbance obtained for VP7 detection must be within the limits established by the Patron Curve.

VII. TEST PROCEDURE

1. All reagents must be allowed to come to room temperature before use.
2. **Sample and standards addition:**
Add to the respective wells 100 µl of negative control and the 5 standards in duplicate. Add 100 µl of, at least, two dilutions of the samples (prepared as it is described above) in duplicate.
3. Seal the plate and incubate **1 hour at 37°C.**
4. Wash 4 times following the described procedure.

5. Add 100 µl of conjugate to each well. Seal the plate and incubate **1 hour at 37°C**.
6. Wash 4 times following the described procedure.
7. Add 100 µl of substrate to each well and incubate for **10 minutes** at room temperature. (Count the time since the first well had been filled). .
8. Add 100 µl of stop solution to each well. ATTENTION: The stop solution must be dispensed in the same order in which the substrate was added.
9. Read the OD of each well at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

Read at 450 nm

A) Validation of the test:

The test will be valid if:

- The absorbance obtained with the point of 50 ng/ml standard is higher than 0.9 **and**
- The absorbance obtained with the negative control is lower than the absorbance obtained with the point of 3.12 ng/ml

B) Interpretation of the results:


The amount of VP7 contained in the sample is quantified in the following way:

- Make the Patron Curve (see the example)
- Determine the concentration of VP7 in the sample extrapolating its OD value in the curve (see the example).

D) Example:


1. Negative control OD = 0.06
2. Standards OD:

ng/ml	Abs450
50	1.38
25	0.65
12,5	0.32
6,25	0.18
3,12	0.11

3. The concentration of VP7 in the sample (C) can be extrapolated in the curve by using the appropriated software. For example:

- ✓ Software: Excel (Microsoft)
- ✓ Select the cells with the absorbance and concentration and clic "Graphic assistance"
- ✓ Select "Type of graphic": "XY (dispersion)" and "finish"
- ✓ Click on "Graphic" and select " Add tendency line"
- ✓ In "Type", select "Lineal"
- ✓ In "Options", select " Show formula in the graphic"
- ✓ In this formula, substitute x for the OD of the sample to obtain "C"

Sample	Dilution(D)	OD	ng/ml (C)*	(C xD)/1000	Mean (µg/ml)
B	1/1000	0,578	21,076	21,08	21.34
	1/2000	0,299	10,80	21,60	

*Obtained by replacing x for the OD of the sample in the tendency line

- Substitute this data (C) in the formula in order to determine the amount of VP7 of the sample.

$$\text{VP7 in } \mu\text{g/ml} = \text{C} \times \text{D} / 1000$$

C: concentration of VP7 contained in the sample determined in the curve

D: is the factor of dilution (100, 500, 1000, etc.)

1/1000: factor from the conversion: ng/ml to µg/ml

C) Additional information

- In order to calculate correctly, the OD values of the samples must be within the limits of the curve. On the contrary, the assay must be repeated with other dilutions of the sample.
- Just in case that the OD value obtained be slightly higher than the point of 50 ng/ml (until 15%) the assay could be considered valid

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in by: