## INGEZIM MAEDI SCREENING

Prod Ref: 13.SLS.K1

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---

## I. TECHNICAL BASIS

Maedi Visna Virus (MVV) is a member of the genus lentivirus of the family Retroviridae. It is a viral disease that affects sheep and goats occasionally. This virus is closely related to the of the caprine arthritis-encephalomyelitis virus (CAEV), a lentivirus, occurring most frequently in goats. Both viruses share many characteristics, and often are considered altogether the lentivirus of small ruminants (SRLV). The SRLV can be classified in A genotype, which includes the Maedi Visna Virus and genotype B assigned to the CAEV. Our kit is capable of detecting specific antibodies of both MVV and CAEV diseases using specific peptides of both genotypes (A and B) and is based on the technique of indirect enzyme

immunoassay, which is described briefly below. On a solid support (polystyrene plate), the preparation of peptides of the genotypes A and B specific MVV and CAEV respectively is adsorbed.

When on the plate the serum sample is dispensed in the case of containing specific antibodies, these will be attached to the plate peptides. After successive washes to remove material not bound, the presence of ovine or caprine immunoglobulins can be detected through the addition of a specific peroxidase-conjugate, which will produce a colorimetric reaction in the presence of the appropriate substrate

---

## II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
10. Substrate and substrate buffer must be handle with care, they are very sensible to light and contamination.

---

## III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

**WARNING:** Stop Solution may precipitate at this temperature, if it happens, warm up at +37°C before use

---

## IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF SAMPLES:

Make the 1/20 dilution of each sample (i.e. 10 µl of serum and 190 µl of diluent)

## VI. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionized water (40 ml of concentrated solution and 960 ml of water). When ready this solution remains stable stored between 2°C and 8°C.

- **Preparation of controls and conjugate:**

Controls and conjugate are ready to use. **Do not dilute.**

- **Preparation of the substrate solution: to make immediately before use.**

It is ready to use. Do not dilute.

## VII. TEST PROCEDURE

1. All reagents must be allowed to bring room temperature before use.
2. Add samples diluted as it's specified in the previous instructions. Add 100 µl of Positive and Negative control serum to their respective wells We recommend to run controls in duplicate wells. **Incubate the plate for 45 min ± 2 min at 37°C**
3. Wash 3 times following the described procedure.
4. Add 100 µl of conjugate to each well and **incubate the plate for 30 min ± 2 min at 37°C.**
5. Wash 3 times following the described procedure.
6. Add 100 µl of substrate to each well. Keep the plate for **15 min at room temperature.**
7. Add 100 µl of stop solution to each well.
8. Read the OD of each well with a spectrophotometer at 405 nm within 5 min after the addition of stop solution.

## VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **405 nm**.

$$S/P = \frac{OD_{405} \text{ Sample} - OD_{405} \text{ C (-)}}{OD_{405} \text{ C (+)} - OD_{405} \text{ C (-)}}$$

### A. Validation criteria:

The test is considered valid when:  
Mean OD of the positive control is > 0.6  
Mean OD of the negative control is < 0.2

### B. Positivity index calculation (S/P):

If you run samples in duplicate, OD values will be calculated for each sample as the arithmetic mean of both values.

To obtain S/P of the sample:

### C. Results Interpretation:

Samples with S/P higher than 0.45 will be considered as Positives.

Samples with S/P lower than 0.40 will be considered as Negatives.

Samples with S/P equal or between both values must be considered as doubtful. We recommend to use the confirmation assay (INGEZIM MAEDI CONFIRMACIÓN) with doubtful and weak positive sera (S/P<0.6)

Desarrollado y fabricado en España por:  
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: ingenasa@ingenasa.com  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in by: