# INGEZIM MAEDI CONFIRMATION

Prod Ref: 13.SLC.K1

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Indirect immunoenzymatic assay for the confirmation of specific detection of serum antibodies to small ruminants Lentivirus (MVV and CAEV) in ovine and caprine samples.

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OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT:

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Distilled or deionised water.  
Micropipettes from 5 to 200 µl.  
Disposable micropipette tips.  
Washing plates device.  
Test tubes from 50 to 250 ml  
ELISA Reader (405 nm filter)

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## I. TECHNICAL BASIS

The Maedi Visna Virus (MVV) is a member of the genus Lentivirus of the Retroviridae family. It is a viral disease that affects sheep and occasionally goats. This virus is closely related to the caprine arthritis-encephalitis virus (CAEV), a Lentivirus, more frequently found in goats. Both viruses share many characteristics, and are often considered altogether as small ruminants lentiviruses (SRLV). The SRLV can be classified in genotype A, which includes the Maedi Visna virus and in genotype B assigned to the CAEV.

This assay has been designed to confirm the positive samples obtained in the screening test (13.SLS.K1). Therefore, in addition to the detection of specific antibodies to both MVV and CAEV through the use of specific peptides of both genotypes (A and B), the kit provides plates coated with negative antigen. The assay is based on an indirect

enzyme immunoassay technique, which is briefly described below.

A preparation of peptides of the genotypes A and B, specific to MVV and CAEV respectively, are adsorbed on a polystyrene plate. On another plate (red plate) the negative peptide is fixed. After adding a serum sample on both plates, in the case of containing specific antibodies, these will be attached to the positive plate. If the antibodies are non-specific they will attach to the positive and to the negative plate. After several washing steps to remove any unbound material, the presence of ovine or caprine immunoglobulins can be detected by the addition of a specific peroxidase-conjugate, which will produce a colorimetric reaction in the presence of an appropriate substrate. The subtraction of the values obtained in both plates will determine whether the antibodies are specific to lentivirus or not.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix reagents or use instructions from different kits.
4. Avoid any contamination of the kit's reagents.
5. Do not use components after expiration dates and do not mix components from different batches.
6. There should be no eating, drinking, or smoking while handling reagents or samples.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Systematically include a positive control and a negative serum each time the assay is run.
10. The substrate and the substrate buffer must be handled with care, as they are very sensitive to light and contamination.

## III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C. WARNING: Stop Solution may precipitate at this temperature, if it happens, warm it up at +37°C before use

## IV. INFORMATION ABOUT THE WASHING PROCEDURE

The washing of the plates can be done using an automatic plate washer, a multichannel pipetting device suitable for dispensing 300 µL on each well, or a squeeze bottle.

After the incubation periods, the washing steps must be done following the instructions below:

- Throw away the content of the plate by briskly turning the plate over, in order to avoid the possible exchange of contents from one well to another.
- Dispense a volume of 300 µl of washing solution in each well.
- Shake the plate gently, avoiding cross-contamination between wells.
- Briskly turn the plate over to empty the wells.
- Repeat the process as many times as indicated in the kit's instructions.
- Prior to emptying the content of the plate in the last washing step, verify that the next reagent to be added is ready to be used. Do not let the plate dry longer than strictly necessary to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use.
- After the last washing step tap the plate upside down on absorbent filter paper to remove any remaining washing solution.

## V. PREPARATION OF SAMPLES:

Dilute the samples at 1/20 dilution (i.e. 15 µl of serum and 285 µl of diluent)

## VI. PREPARATION OF REAGENTS

- **Washing solution:**  
Dilute one part of the concentrated washing solution provided in the Kit with 24 parts of distilled or deionized water (40 ml of concentrated solution and 960 ml of water). Once prepared this solution remains stable when stored between at +.
- **Preparation of the substrate solution**  
Substrate is ready to use. Do not dilute.
- **Preparation of controls and conjugate:**  
Do not dilute. Controls and conjugate are ready to use.

## VII. TEST PROCEDURE

1. Prior to starting the test, bring all reagents to room temperature (22-25°C).
2. **IMPORTANT: All samples must be assayed in the positive and the negative wells.** Add the samples diluted as specified in the previous instructions. Then add 100 µl of the Positive and Negative control serum to their respective wells. We recommend (no mandatory) running controls in duplicate in both wells. **Incubate the plate for 45 min ± 2 min at 37°C.**
3. Wash 3 times following the procedure previously described.
4. Add 100 µl of conjugate to each well and **incubate the plate for 30 min ± 2 min at 37°C.**
5. Wash 3 times following the procedure previously described.
6. Add 100 µl of the substrate to each well. Keep the plate for **15 min at room temperature.**
7. Add 100 µl of the stop solution to each well.
8. Read the absorbance of each well with a spectrophotometer at 405 nm within the following 5 min after the addition of the stop solution.

## VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **405 nm**.

The values obtained in the negative wells must be subtracted from the values obtained in the positive wells

### A. Validation criteria:

The test is considered valid when:

The OD of the positive control in the positive well minus the OD in the negative well is  $\geq 0.6$

The OD of the negative control in the positive well minus the OD in the negative well is  $\leq 0.2$

### B. Positivity index calculation

To obtain the index of the sample:

$$(\text{OD Sample of pos. well} - \text{OD sample of neg. well}) \times 100$$

### C. Results Interpretation:

Samples with an index  $\geq 25$  must be considered as Positive to SRLV antibodies.

Samples with an index  $< 25$  will be considered as Negative to SRLV antibodies.

Desarrollado y fabricado en España por:  
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: ingenasa@ingenasa.com  
[www.ingenasa.es](http://www.ingenasa.es)

IT-73840  
IT-73780

ISO 14001:2015  
9191.INGE

ISO 9001:2015  
9175.ING2

Distributed in

by: