INGEZIM ROTA BOVINO

Prod Ref: 12.RT.K1

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I. TECHNICAL BASIS

The kit is based on an indirect enzymatic immunoassay (**Indirect ELISA**). We make a brief description of the technique below:

We fix the antigen on a solid support (Polystyrene plate). When a sample serum contains specific antibodies against the virus, they will bind to the antigen absorbed on plate. After washing to eliminate all non-fixed material from the sera sample, we can detect presence of bovine immunoglobulins using a

specific Monoclonal antibody conjugated with peroxidase. After addition of the substrate, a colorimetric reaction will appear which could be measured by a spectrophotometer.

In this way, the presence of colour means the presence of antibodies against the virus in the bovine sera and the absence of colour, the absence of specific antibodies.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the kit's reagents.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
10. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
11. Substrate must be handled with care, it is very sensible to light and contamination.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

Once opened, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

The sera samples need to be at 1/100 dilution in diluent to be tested (5 µl of serum in 0.5 ml of diluent). For titration, we recommend to assay 8 two fold

dilutions for each sample (from 1/100 to 1/12800). For your facilities this step could be done directly on the plate

VI. PREPARATION OF REAGENTS

• **Washing solution:**

Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionized water (40 ml of concentrated solution and 960 ml of water). When ready this solution remains stable stored at 4°C.

• **Preparation of control sera (+) and (-):**

Make dilution 1/100 (as sample sera) in diluent (5 µl of serum in 0,5 ml of diluent). If you need to do a partial use of the kit, we recommend to make aliquots of the control sera and store them at -20°C till use.

• **Preparation of Diluent:**

Dilute one volume of the concentrated diluent provided in the kit with 4 volumes of distilled or deionized water. Store this solution at 4°C.

• **Preparation of the conjugate: to make immediately before use.**

Dilute 1/100 with diluent 1x prepared as it's described at point "Preparation of diluent":

- The necessary and sufficient quantity of conjugate for a complete plate is 110 µl of conjugate in 11 ml of diluent.
- The necessary and sufficient quantity of conjugate for an 8 wells strip is 10 µl of conjugate in 1 ml of diluent.

VII. TEST PROCEDURE

1. All reagents (except conjugate) must be allowed to come to room temperature before use.
2. Add 100 µl of each dilution of positive control serum, diluted as it's specified in the previous instructions, to each well of the first row of the plate and 100 µl of each dilution of the negative control serum to the wells of the second row of the plate. Add 100µl of each dilution of sera samples to test, prepared following the previous instructions, on the remainder wells of the plate (use one row per serum). Seal the plate and incubate for **60 min at 37°C**.

This assay can be used as **screening assay**. In that case, controls, as well as samples sera must be tested in two wells at a single dilution of 1/100.
3. Wash 3 times following the procedure previously described.
4. Add 100 µl of conjugate, prepared following the previous instructions, to each well. Seal the plate and incubate **for 30 min at 37°C**.
5. Wash 5 times following the described procedure.
6. Add 100 µl of substrate. Keep the plate for **10 min at room temperature** in darkness.
7. Add 100 µl of stop solution to each well.
8. Read the absorbance of each well with a spectrophotometer at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm**.

A. Validation criteria:

The assay will be valid only when:

- OD of positive control is higher than 0.8
- OD of negative control is lower than 0.2

B. Cut off calculation:

Cut off = OD of negative control plus 0,4

C. Results Interpretation:

• **For screening proposes:**

The OD value must be the average between the OD values obtained for each sample.

- All samples with an OD value higher than Cut off must be considered positive.
- All samples with an OD value lower than Cut off, must be considered negative.

• **For titration purposes:**

The titre of each sample, will be the highest dilution which develops a positive result (OD value higher that Cut off value)

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by:

Empty box for distributor information.