## INGEZIM IBR Compac 2.0

Prod Ref: 12.BHV.K3

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### 12.BHV.K3

## I. TECHNICAL BASIS

The Kit has been designed to detect specific antibodies to BoHV1 in cattle sera and milk samples. This kit is based on a blocking enzymatic immunoassay (Blocking Elisa). The technique is briefly described below:

The BoHV-1 antigen is fixed on a solid support (polystyrene plate). The samples are added to each well. After an incubation period, a specific monoclonal antibody (peroxidase conjugated) specific to gB protein of the BoHV-1 is added. If the serum sample does not contain antibodies against the virus, the conjugate will bind to the plate, whereas if it contains specific antibodies to BoHV-1, these

antibodies will bind to the antigen on the plate, blocking the binding of the conjugate, which will not occur.

After washing the plate to eliminate all non-fixed material, the presence or absence of labelled conjugate can be detected by adding a substrate which, in presence of the peroxidase, will develop a colorimetric reaction.

Due to the high specificity of the monoclonal used as the conjugate, the assay has high specificity and sensitivity indexes.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix reagents or use instructions from different kits.
4. Avoid any contamination of the kit reagents.
5. Do not use components after expiry dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking while handling reagents or samples.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Systematically include a positive control and a negative serum each time the assay is run.
10. The stop solution is a strong acid and therefore must be handled with care. In case of an accidental contact with skin, rinse gently with water.
11. The substrate must be handled with care, as it is very sensitive to light and contamination.

## III. STORAGE OF COMPONENTS

All reagents and components must be stored between 2°C and 8°C.

**Once opened**, the control sera are stable for one month between +2°C and +8°C. If they are not going to be used during this period, we recommend storing them below -20°C

## IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps can be done using an automatic plate washer or a multichannel pipetting device suitable for dispensing 300 µl in each well.

After the incubation periods, the washing steps must be done following the instructions below:

- Throw away the content of the plate by briskly turning the plate over, in order to avoid the possible exchange of contents from one well to another.
- Dispense a volume of 300 µl of washing solution in each well.
- Shake the plate gently, avoiding cross-contamination between wells.
- Briskly turn the plate over to empty the wells.
- Repeat the process as many times as indicated in the kit instructions.
- Prior to emptying the content of the plate in the last washing step, verify that the next reagent to be added is ready to be used. Do not let the plate dry longer than strictly necessary.
- After the last washing step tap the plate upside down on absorbent filter paper to remove any remaining washing solution.

## V. PREPARATION OF SAMPLES:

### **Sera:**

Samples must be diluted 1/2 in sample diluent DE01-01. The dilution can be done in the same well by first putting 50 µl of diluent and secondly 50 µl of serum .

### **Milk (individual and tank)**

Must be tested undiluted. Fresh, refrigerated or previously frozen milk may be tested. To eliminate interference from lipids, partially skimmed milk should be either by centrifugation (15 min at 2000xg) or leaving overnight at 4°C to observe the formation of a lipid layer on the surface. Sample should be collected under the layer of lipids and that will be the one used in the assay. To obtain whey, freeze/thaw the milk sample

## VI. PREPARATION OF REAGENTS

• **Washing solution:**

Dilute one part of the concentrated washing solution provided in the Kit with 24 parts of distilled or deionized water (40 ml de concentrated solution and 960 ml of water). Once prepared, this solution remains stable when stored between +2°C and +8°C.

• **Preparation of controls:**

Controls must be diluted ½ in diluent DE01-1 if the samples are sera, and diluted 1/5 if the samples are milk.

ATTENTION! Once opened, the control sera are stable for one month when stored between +2°C and +8°C. For long term storage, we recommend making aliquots and storing them at -20°C.

• **Preparation of the conjugate: (to be prepared immediately before use):**

Dilute 1/100 the required volume of conjugate with the diluent provided with the kit:

The necessary and sufficient amount of conjugate for a complete plate is 110 µl of conjugate with 11 ml of diluent. The necessary and sufficient amount of conjugate for an eight wells strip is 10 µl of conjugate with 1 ml of diluent.

Shake the solution well before use. Prepare only the amount of conjugate necessary to run the test, as the remaining solution must be discarded.

## VII. TEST PROCEDURE

1. Prior to starting the test, bring all reagents to room temperature (20-25°C) (except conjugate).
1. **Adding the samples (serum or milk):**  
Dispense 100 µl/well of each sample prepared as indicated in section V. Add last dilute controls as indicated in section VI. It's recommended to evaluate samples and controls in Seal the plate and **incubate for 1 hour at 37°C if the samples are sera and 16 - 20 h at room temperature in humidity chamber if the samples are milk**
2. Wash 5 times following the procedure previously described.
3. Add 100 µl of the conjugate (prepared following previous instructions) to each well. Shake the plate carefully to help the correct homogenisation of the components mixed on each well. Seal the plate and **incubate for 30 min at +37°C.**
4. Wash 5 times following the procedure previously described.
6. Add 100 µl of the substrate solution, to each well. Keep the plate for **15 min at room temperature, in a dark place.**
7. Add 100 µl/well of the stop solution.
8. Read the OD of each well at 450 nm within the following 5 min after the addition of the stop solution.

## VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm.**

### A. Validation criteria:

- OD of Negative control >0.750
- To validate the assay, the ratio:

OD (positive Control) / OD (negative control)  
must be lower than 0.3

### B. Results Interpretation:

When running the test in duplicate, the OD of the sample will be calculated as the arithmetic mean of the OD values in both wells.

The blocking percentage for each sample must be calculated as follows:

$$\text{blocking \%} = \frac{\text{OD neg Control} - \text{OD sample}}{\text{OD neg Control} - \text{OD pos Control}} \times 100$$

### Serum

- Samples with a blocking % lower than 25% must be considered as Negative.
- Samples with a blocking % higher than 30% will be considered as Positive.
- Samples with blocking % between both should be considered as doubtful. We recommend re-testing after 3 weeks.

### Milk

- Samples with a blocking % lower than 40% must be considered as Negative.
- Samples with a blocking % higher or equal than 40% will be considered as Positive.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in by:

IT-73840  
IT-73780

ISO 14001:2015  
9191.INGE

ISO 9001:2015  
9175.ING2