## INGEZIM CIRCO IgG

Prod Ref: 11.PCV.K1

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## 11.PCV.K1

### I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Elisa). We make a brief description of the technique below:

The antigen (PCV recombinant protein) is fixed in a solid support (polystyrene plate). When a sample serum contains specific antibodies against the virus, they will bind to the antigen absorbed on plate while if the serum sample does not contain specific antibodies they will not bind the antigen. If we add a specific monoclonal antibody against the porcine immunoglobulins (conjugated with peroxidase), it will

bind the antibodies of the serum. If the serum sample contains specific antibodies, conjugate will bind them, whereas if it does not contain specific antibodies the conjugate will be eliminated in the following washing step. After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled conjugate by adding the specific substrate that in presence of the peroxidase will develop a colorimetric reaction.

The antigen coated to the plate in our kit consists of a soluble recombinant protein from the virus.

### II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions for use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiry dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Substrate and substrate buffer must be handled with care.
10. For each utilisation of the Kit use a fresh substrate preparation.
11. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.

### III. STORAGE OF COMPONENTS

All reagents and components must be stored entre +2°C y +8°C.

### IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

### V. PREPARATION OF SAMPLE

Samples to be used in the assay must be porcine sera.

1/200 dilution of the sera must be done in diluent (5 µl of sera and 1 ml of diluent).

### VI. PREPARATION OF REAGENTS

#### • **Washing solution:**

Dilute one part of the concentrate washing solution provided in the kit with 24 parts of distilled or deionised water (i.e. mix 40 ml of concentrate with 960 ml of water to prepare 1L of working washing solution). When ready this solution remains stable at +4°

#### • **Control sera:**

Controls are supplied ready to use. Do not dilute

#### • **Diluent:**

Dilute one part of the concentrate diluent with 4 parts of distilled or deionised water (mix 100 ml of

concentrate with 400 ml of water). When ready this solution remains stable at +4°

• **Conjugate:**

Conjugate is supplied ready to use. Do not dilute

## VII. TEST PROCEDURE

- All reagents must be allowed to come to room temperature before use.
- Add 100 µl of positive control serum to two wells (e.g. A1 and B1), and 100 µl of the negative control serum to other two wells (e.g. A2 and B2). Add 100 µl of sera samples to test (prepared as specified in point VI), on each remainder wells. We recommend the use of two wells per sample. **Seal the plate and incubate for 1 h at room temperature.**
- Wash four times as specified in point V.
- Add 100 µl of conjugate to each well Cover the plate and **incubate the plate for 30 minutes at room temperature.**
- Wash 6 times following the described procedure.
- Add 100 µl of substrate, to each well. Keep the plate **for 10 min at room temperature.**
- Add 100 µl of stop solution, to each well.
- Read the OD of each well at 450nm within 5 min after the addition of stop solution.

## VIII. READING AND INTERPRETATION OF THE RESULTS

### 1. Validation of the test

The test could be considered valid when:

- ⇒ The OD of the **Negative Control** (NC) is lower than **0,35**
- ⇒ The OD of the **Positive Control** (PC) is higher than **0,7**

### 2. Interpretation:

If samples are analyzed in duplicate, OD of the sample will be calculated as the arithmetic mean of OD values in both wells.

• **Cut Off calculation:**

The cut off must be calculated as follow:

- Negative Cut off = OD negative control+0.2
- Positive Cut off = OD negative control+0.25

Samples with OD higher than positive cut off value should be considered positive to PCV antibodies. Samples with OD lower than negative cut off value should be considered negative. OD between both cut off will be considered doubtful.

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The titre of positive samples will be determined as follow:

- a) Calculate the S/P

S/P = OD of sample / OD of positive control

- b) Titre is calculated as follow

$$\text{Titre} = 53 (e^{3.2 \times x})$$

Where, **e** is the natural logarithm base and **x** is the S/P of the sample

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INMUNOLOGIA Y GENETICA APLICADA, S.A.