## INGEZIM Circovirus IgG/IgM

Prod Ref: 11.PCV.K2

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## 11.PCV.K2

### I. TECHNICAL BASIS

This kit is based on a capture enzymatic immunoassay. We make a brief description of the technique below:

Two monoclonal antibodies (Mabs) specific for pig IgM and IgG are fixed on different plates. When a pig serum is added on each plate, the IgMs and IgGs present in the sample are captured by the Mab adsorbed on the plate. After washing to eliminate all non fixed material from the serum sample, we add the viral antigen (recombinant empty capsid) which is captured by the IgM and/or IgG only if the pig serum sample contained specific immunoglobulin M or G to PCV.

After additional washing we add a peroxidase-conjugated Mab specific to PCV2. After washing and addition of the substrate a colorimetric reaction will be developed which could be measured by an ELISA reader.

In this way the presence of colour means the presence of antibodies IgM or IgG against the virus in the pig sera, and the absence of colour the absence of specific antibodies.

The aim of this kit is to provide the users with a reliable diagnostic technique for this disease.

### II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use. **¡IMPORTANT!**
3. Do not mix instructions or reagents from different kits.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
6. Avoid any contamination of the reagents
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Substrate is very sensitive to light and contamination. Stop solution is a strong acid. Handle with care.
10. For each utilisation of the kit, control positive and negative sera must be tested in a systematic way.

### III. STORAGE OF COMPONENTS

All reagents and components must be stored between +2°C and +8°C.

### IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

### V. PREPARATION OF SAMPLE

The sera samples needs to be at **1/100 dilution** in diluent to be tested.

### VI. PREPARATION OF REAGENTS

#### • **Washing solution:**

Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionized water (i.e. 100ml of

concentrate and 2,4 L of water). When ready, this solution remains stable between +2°C and +8°C.

- **Preparation of the control sera (+) and (-) :**  
Controls are ready to use. Do not dilute. There are a positive control for IgM, a positive control for IgG and only a negative control for both.

- **Preparation of the recombinant viral antigen:**

Antigen is ready to use. Do not dilute

- **Preparation of the conjugate : to make immediately before use.**

Dilute 1/100 with diluent.

- The necessary and sufficient quantity of conjugate for a complete plate is: 110 µl of conjugate in 11 ml of diluent.
- The necessary and sufficient quantity of conjugate for a 8 wells strip is: 10µl of conjugate in 1 ml of diluent.

Prepare only the quantity needed for each time because the remainder volume has to be rejected.

Shake very well the solution prior to use.

## VII. TEST PROCEDURE

1. All reagents and wells to be used must be allowed to come to room temperature before use (except conjugate).
2. Add 100 µl of diluted sera samples and positive and negative controls to IgM and IgG wells (the negative control to both and the positive controls to corresponding wells). We recommend to run the test (samples sera and controls), in duplicates. Seal the plate and **incubate for 60 minutes at 37°C**.
3. Wash 4 times following the procedure previously described.
4. Add 100 µl of antigen to each well. Seal the plate and **incubate 30 minutes at 37°C**
5. Wash 4 times following the described procedure
6. Add 100 µl of conjugate, prepared following the previous instructions, to each well. Seal the plate and **incubate for 30 minutes at 37°C**.
7. Wash 4 times following the described procedure.
8. Add 100 µl of substrate solution, to each well. Keep the plate for **5 min at room temperature**.
9. Add 100 µl of stop solution to each well.
10. Read the absorbances with an ELISA reader at 450 nm within 5 min after the addition of stop solution.

## VIII. READING AND INTERPRETATION OF THE RESULTS

The reading must be done at **450 nm**.

### 1. Validation of the test

The test is valid when:

- **O.D positive control in IgG and IgM wells are higher than 0.7**
- **O.D negative control in IgG and IgM wells are lower than 0,3**

### 2. Results Interpretation:

Using the corresponding positive control value, the cut off points for IgG and IgM must be calculated as follows:

- **CUT OFF: IgM**  
O.D 450 nm IgM positive control x 0.4
- **CUT OFF: IgG**  
O.D 450 nm IgG positive control x 0.3

If samples in duplicate wells are run, the OD value for the sample will be the mean of both values.

- **SAMPLES WITH OD VALUES < CUT OFFS**  
These kind of samples have titre lower or equal than 1/100, (for IgG and/or IgM). These values corresponds with negative animals or animals with a very old contact with virus
- **SAMPLES WITH OD VALUES > CUT OFFS:**  
In this case you could distinguish three kind of samples:

- IgM > IgG OD value: Active infection
- IgG > IgM OD value and IgM titre > 1/100: Recent infection (between one or two months)
- IgG > IgM OD value and IgM titre < 1/100: Old infection by Circovirus.

Desarrollado y fabricado en España por:  
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tif: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in \_\_\_\_\_ by: