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título =  $40(e^{1,9 * 1,22}) = 406$

## I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Elisa). We make a brief description of the technique below:

The antigen (ADV proteins) is fixed in a solid support (polystyrene plate). When a sample serum contains specific antibodies against the virus, they will bind to the antigen absorbed on plate while if the serum sample does not contain specific antibodies they will not bind the antigen. If we add a specific monoclonal antibody against the porcine immunoglobulins (conjugated with peroxidase), it will bind the antibodies of the serum. If the serum sample

contains specific antibodies, conjugate will bind them, whereas if it does not contains specific antibodies the conjugate will be eliminated in the following washing step. After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled conjugate by adding the specific substrate that in presence of the peroxidase will develop a colorimetric reaction.

The antigen coated to the plate in our kit consists of a soluble protein extract from the virus containing all the antigenic proteins present in vaccine strains as well as in field strain.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or Smoking where specimens or Kit reagents are being handle
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Substrate must be handled with care because it is very sensitive to light and contamination
10. Stop solution is a strong acid. Handle with care.
11. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.

## III. STORAGE OF COMPONENTS

Store all plates and reagents between +2°C and +8°C.

## IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF SAMPLES:

Samples to be used in the assay must be porcine sera.  
1/25 dilution of the sera must be done in diluent (10 µl of sera and 240µl of diluent).

## VI. PREPARATION OF REAGENTS

### • **Washing solution:**

Dilute one part of the concentrate washing solution provided in the kit with 24 parts of distilled or deionized water. When ready this solution remains stable at +4°

### • **Control sera:**

Controls are supplied ready to use. Do not dilute

### • **Diluent:**

Diluent is supplied ready to use. Do not dilute.

### • **Preparation of the conjugate: to make immediately before use.**

Dilute the needed quantity of conjugate provided in the kit 1/100 into diluent DE14-01 (e.g.110 µl of conjugate into 11ml of diluent is enough quantity for one plate).

Shake very well the solution before the use. Prepare only the quantity needed for each time because the remainder volume has to be rejected.

## VII. TEST PROCEDURE:

- All reagents (except conjugate) must bring to room temperature before use.
- Add 100 µl of positive control serum to two wells (e.g. A1 and B1), and 100 µl of the negative control serum to other two wells (e.g. A2 and B2). Add 100 µl of sera samples to test (prepared as specified in point V), on each remainder wells. We recommend the use of two wells per sample. Seal the plate and incubate for **1 h at 37°C**.
- Wash four times as specified in point V.
- Add 100 µl of specific conjugate (prepared following previous instructions) to each well. Cover the plate and incubate the plate for **30 minutes at 37°C**.
- Wash 5 times following the described procedure.
- Add 100 µl of substrate, to each well. Keep the plate for **5 min at room temperature**.
- Add 100 µl of stop solution, to each well.
- Read the OD of each well at **450nm** within 5 min after the addition of stop solution.

## VIII. READING AND INTERPRETATION OF THE RESULTS:

### 1. Validation of the test

The test could be considered valid when:

- The OD of the Negative Control (NC) is lower than 0.2
- The OD of the Positive Control (PC) is higher than 0,5

### 2. S/P Calculation:

When you are running duplicates, OD of the sample will be calculated as the arithmetic mean of OD values in both wells.

- a) Calculate the S/P

$$S/P = \frac{S}{PC}$$

Where,  
S is the sample OD  
PC is the positive control OD

### 3. Interpretation :

- Samples with s/p higher Than 0.35 should be considered positive-
- Samples with s/p lower than 0.35 should be considered negative.

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- The titre of positive samples will be determined as follow:

$$\text{Titre} = 40 * (e^{1.9x})$$

Where, **e** is the natural logarithm base and **x** is the S/P of the sample.

---

## 5. Example:

---

	<u>DO</u>
<b>Control +</b>	<b>0,743</b>
<b>Control -</b>	<b>0,102</b>
Sample 1	0.910

### 1. Validation:

Positive control OD value = 0,743 (> 0,5)

Negative control OD value = 0,102 (< 0,2)

### 2. S/P Calculation

$$s/p \text{ Sample} = 0.910/0.743 = 1,22$$

### 3. Interpretation

Positive Sample ( > 0,35)

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a) S/P (x)

$$x = 0,910/ 0,743 = 1.22$$

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$$\text{Titre} = 40 (e^{1.9 \times 1,22}) = 406$$

Desarrollado y fabricado en España por:  
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tif: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in \_\_\_\_\_ by: