INGEZIM ADV gE PLUS

Prod Ref: 11.GEP.K3

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I. TECHNICAL BASIS:

The gE test uses a monoclonal antibody (Mab) specific for gE glycoprotein of ADV and microtiter plates coated with ADV antigen (gE). After sample addition, if antibodies are absent in the test serum, the gE antigen will remain free and the Mab peroxidase-conjugated specific for gE, will bind. After substrate

addition a colour reaction is produced. When gE antibodies are present in the test sample they will bind to the gE antigen of the well and will block the binding of the conjugated Mab. In this case no colour reaction will be produced.

II. PRECAUTIONS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Avoid any contamination of the reagents of the Kit.
4. Do not use components after expiration date.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handle.
6. Do not pipette by mouth.
7. Use a new tip for each serum sample.
8. Substrate must be handle with care, it is very sensible to light and contamination
9. Stop solution is a strong acid. Handle with care
10. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.

III. STORAGE OF COMPONENTS:

All reagents and components must be stored between 2°C and 8 ° C.

Once opened, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

IV. INSTRUCTIONS FOR WASHING STEPS:

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well. After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate.
- Dispense a volume of 300µl of washing solution on each well.

- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

The sera sample needs to be at ½ dilution to be tested. The dilution of sera samples could be done directly on the wells by adding 50 µl of serum and 50 µl of sera diluent to each well.

VI. REAGENTS PREPARATION INSTRUCTIONS:

• Washing solution:

Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionized water (40 ml of concentrate washing solution into 960 ml of distilled water). This solution remains stable when is stored between +2°C and +8°C.

• Positive and Negative Controls:

Controls need to be at 1/2 dilution to be tested. The dilution of control samples could be made directly on the wells by adding 50 µl of serum and 50 µl of sera diluent to each well.

VII. TEST PROCEDURE:

1. All reagents must be allowed to come to room temperature before use.
2. Add 50 µl of diluent and 50 µl serum (controls and samples) to each well. Shake carefully to homogenate the mixture, seal the plate and incubate for 2 hours at R.T. (20-25°C) (preferably with agitation) or overnight between +2°C and +8°C.
3. Wash 3 times following the described procedure.
4. Add 100 µl of conjugate to each well. Seal the plate and incubate for 30 minutes at R.T. (20-25 °C)
4. Wash 4 times following the procedure.
5. Add 100 µl of substrate, to each well. Keep the plate for 10-15 min at room temperature in a dark place.
6. Add 100 µl of stop solution to each well.
7. Read the OD at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION:

1. TEST VALIDATION

- P/N ratio must be lower than 0,3:
P= OD of positive control
N= OD of negative control
- Negative control's OD must be higher than 0,800.

2. CUT OFF CALCULATION:

Positive CUT OFF= $N - [(N-P) \times 0,4]$
Negative CUT OFF= $N - [(N-P) \times 0,35]$

N= OD of Negative control (average)
P= OD of positive control (average).

3. RESULTS INTERPRETATION:

If you are running samples in duplicate wells, Optical Density (OD) value must be calculated as the arithmetic average of the two wells. Sample is positive when its OD value is lower than or equal to the Positive cut off value. Sample is negative when its OD value is higher than or equal to the negative cut off value. Samples with OD values between both cut off values must be scored as doubtful. In these cases we recommend to run a second sample of the same animal collected 2 weeks later.

4. INHIBITION PERCENTAJE CALCULATION (IP_M)

$$IP_M = \frac{OD_{CN} - OD_{MUESTRA}}{OD_{CN} - OD_{CP}} \times 100$$

IP_M ≥ 40: Sample positive to gE antibodies

IP_M ≤ 35: Sample negative to gE antibodies

40 ≥ IP_M ≥ 35: Doubtful sample to gE antibodies

5. PRACTICAL EXAMPLE:

- **Obtained OD average values:**

NC = 1.160 (N); PC = 0.110 (P)

Sample 1 = 0.373; Sample 2 = 1.395

- **Assay validation:**

P/N ratio = 0,110 / 1.160 < 0.3

1,160 > 0,800

- **CUT OFF Calculation:**

Cut Off (+) = $1,160 - \{(1,160 - 0,110) \times 0,4\} = 0,740$

CUT OFF (-) = $1,160 - \{(1,160 - 0,110) \times 0,35\} = 0,792$

- **Results:**

Sample 1 = Positive for gE

Sample 2 = Negative for gE

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: