Desarrollado y fabricado en España por:  
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tif: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: ingenasa@ingenasa.com  
www.ingenasa.es

Distributed in \_\_\_\_\_ by:

## INGEZIM INFLUENZA PORCINA

Prod Ref: 11.FLU.K1

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Enzymatic immunoassay for the  
detection and quantification of  
specific antibodies for influenza A  
viruses in porcine serum

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Distilled or deionised water.  
Micropipettes from 5 to 200 µl.  
Disposable micropipette tips.  
Washing plates device.  
Test tubes from 50 to 250 ml  
ELISA Reader (405 nm filter)

**Preparation of the conjugate (to be prepared immediately before its use):**

Dilute 1/100 the required volume of conjugate with the diluent provided in the kit (DE02):

→The necessary and sufficient amount of conjugate for a complete plate is: 110 µl of conjugate with 10.9 ml of diluent.

→The necessary and sufficient amount of conjugate for an 8 well strip is: 10µl of conjugate with 990 µl of diluent.

*Shake the solution well before use. Dilute only the volume of conjugate necessary to run the test, as the remaining solution must be discarded.*

**VI. PREPARATION OF SAMPLES:**

Sera samples must be assayed at 1/100 dilution (i.e. 5µl of serum with 495 µl of diluent)

**VII. TEST PROCEDURE**

- Prior to starting the test, bring all reagents to room temperature (22-25°C) (except conjugate).
- Add the samples, diluted as specified in the previous instructions, to the wells of the plate. Add 100 µl of the positive and negative control serum to their respective wells. We recommend running samples and control in duplicate. **Incubate the plate for 1 hour at 37°C.**
- Wash 4 times following the procedure previously described.
- Add 100 µl of the conjugate (prepared according to the previous instructions) to each

well. **Incubate the plate for 1 hour at room temperature (+18-25°C).**

- Wash 5 times following the procedure previously described.
- Add 100 µl of the substrate to each well. Keep the plate in darkness for **10 min at room temperature.**
- Add 100 µl of the stop solution to each well.
- Read the OD of each well at 450 nm within the following 5 min after the addition of the stop solution.

**VIII. READING AND RESULT INTERPRETATION**

The reading must be done with a spectrophotometer at **450 nm.**

**A. Validation of the test:**

The test is considered valid when:

- The ratio:  
$$\frac{\text{OD Negative Control}}{\text{OD Positive Control}} < 0.2$$
- OD of positive control  $\geq 1$

**B. Interpretation of the results:**

Calculate the ratio: S/P

Sample OD

Positive control (OD)

- Samples with S/P values  $\geq 0,2$  are positive to antibodies to influenza A viruses.
- Samples with S/P values  $< 0,2$  are considered negative to antibodies to Influenza A viruses.

**C. Titer of the samples:**

The titre of a sample must be calculate as follows:

$$T(\text{Titter}) = 2222 \times [(S/P)^{1,93}]$$

**D. Example:**

OD Negative control = 0,162  
OD Positive control = 1,111  
OD Sample = 0,921

Validation =  $0,162 / 1,111 = 0,145 (< 0,2)$   
SP ratio =  $0,921 / 1,111 = 0,829$   
Sample titre =  $2222 \times (0,829^{1,93}) = 1547$

## I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). The technique is briefly described below:

The antigen is fixed on polystyrene plates. When a sample serum contains specific antibodies against the virus, they will bind to the antigen adsorbed on the plate. After a washing step to eliminate all non-fixed material the presence of swine antibodies can be

detected using a specific peroxidase conjugate. After the addition of the substrate a colorimetric reaction will appear which can be measured by a spectrophotometer.

The development of colour indicates the presence of antibodies against the virus in the swine sera, whilst the absence of colour absence of colour means that the sample is negative to specific antibodies.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix reagents or use instructions from different kits.
4. Avoid any contamination of the kit reagents.
5. Do not use components after expiry dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking while handling reagents or samples.
7. Do not pipette by mouth.
7. Use a new tip for each serum sample.
8. The substrate must be handled with care, as it is very sensitive to light and contamination.
9. Systematically include a positive control and a negative serum each time the assay is run.
11. The stop solution is a strong acid and therefore must be handled with care.

## III. STORAGE OF COMPONENTS

Plates and reagents must be stored between +2°C y +8°C.

## IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps can be done using an automatic plate washer or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the instructions below:

- Throw away the content of the plate by briskly turning the plate over, in order to avoid the possible exchange of contents from one well to another.
- Dispense a volume of 300 µl of washing solution in each well.

- Shake the plate gently, avoiding cross-contamination between wells.
- Briskly turn the plate over to empty the wells.
- Repeat the process as many times as indicated in the kit instructions.
- Prior to emptying the content of the plate in the last washing step, verify that the next reagent to be added is ready to be used. Do not let the plate dry longer than strictly necessary.
- After the last washing step tap the plate upside down on absorbent filter paper to remove any remaining washing solution.

## V. PREPARATION OF REAGENTS

### • **Washing solution:**

Dilute one part of the concentrated washing solution provided in the kit with 24 parts of distilled or deionized water (40 ml of concentrated solution and 960 ml of water). Once ready this solution remains stable at +4°C.

### • **Control Sera:**

Control sera are ready to use. Do not dilute

### **Diluent:**

Dilute one part of the concentrated diluent with 4 parts of distilled or deionized water (100 ml of concentrated and 400 ml of water). Once prepared, this solution remains stable when stored between +2°C and +8°C

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