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## I. TECHNICAL BASIS

The kit is based on a Double antibody sandwich enzymatic immunoassay (DAS or capture Elisa). This technique has great advantages over other techniques which have been used traditionally: higher sensitivity and specificity rates, faster and easily automatable. The technique is briefly described below:

A specific monoclonal antibody (Mab) against PPV is fixed on a solid support (polystyrene plate). When a sample containing viral antigen is added, the viral particles will bond to the Mab. After washing the plate to remove all non-fixed material, we add a specific conjugate: a biotinylated monoclonal antibody to PPV.

In presence of the viral antigen this Mab will bind to it. After a second washing step, we add a second conjugate (Streptavidin-peroxidase) that will bind only to the biotin present in each well in the case the antigen had been previously bound. We can demonstrate the presence of this second conjugate by adding a chromogen substrate which, in presence of the peroxidase, will produce a colorimetric reaction.

The use of these very specific and stable Mabs guarantees the objectivity, specificity, safety and accuracy of the assay.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix reagents or use instructions from different kits.
4. Avoid any contamination of the kit reagents.
5. Do not use components after expiry dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking while handling reagents or samples.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. The substrate and the substrate buffer must be handled with care.
10. Systematically include a positive control and a negative serum each time the assay is run.
11. The stop solution may precipitate between +2°C and +8°C. When it happens, warm it up at 37°C before use.

## III. STORAGE OF COMPONENTS

Store all plates and reagents between +2°C and +8°C.

## IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps can be done using an automatic plate washer or a multichannel pipetting device suitable for dispensing 300 µl in each well.

After the incubation periods, the washing steps must be done following the instructions below:

- Throw away the content of the plate by briskly turning the plate over, in order to avoid the possible exchange of contents from one well to another.
- Dispense a volume of 300 µl of washing solution in each well.
- Shake the plate gently, avoiding cross-contamination between wells.
- Briskly turn the plate over to empty the wells.

- Repeat the process as many times as indicated in the kit instructions.
- Prior to emptying the content of the plate in the last washing step, verify that the next reagent to be added is ready to be used. Do not let the plate dry longer than strictly necessary.
- After the last washing step tap the plate upside down on absorbent filter paper to remove any remaining washing solution.

### **ATTENTION!**

**You are working with biological material that could contain infectious agents. All rejected material needs to be sterilised. We recommend using a recipient containing NaOH to throw away the liquids after washing the plate.**

## V. PREPARATION OF REAGENTS

---

- **Washing solution:**

Dilute one part of the 25x concentrated washing solution provided in the Kit with 24 parts of distilled or deionized water. Once prepared, this solution remains stable when stored between +2°C and +8°C.

- **Preparation of the conjugates I and II: to be prepared immediately before use.**

For each conjugate do the preparation process separately. Dilute the required volume of conjugate 1/100 with the diluent provided in the kit:

➤ The necessary and sufficient amount of conjugate for a complete plate is 110 µl of conjugate with 11 ml of diluent.

➤ The necessary and sufficient amount of conjugate for an eight wells strip is 10 µl of conjugate with 1 ml of diluent.

Shake the solution well before use.

Prepare only the amount of conjugate necessary to run the test, as the remaining solution must be discarded.

- **Positive and Negative controls:**

Prepare a 1/100 dilution of the controls with the diluent provided in the kit. (5 µl of control with 500 µl with diluent).

- **Substrate: It is supplied ready to use**

## VI. PREPARATION OF SAMPLES:

---

Macerate 1 gr. of tissue sample in 10 ml of diluent. After a good homogenisation of the sample, centrifuge the homogenate at 1500 r.p.m. for 10 min.

The supernatant obtained will be used as the assay sample.

## VII. TEST PROCEDURE

---

1. Prior to starting the test, bring all reagents to room temperature (22-25°C) (except conjugate).
2. Add 100 µl of the positive control and 100 µl of the negative control (diluted as specified above) to two wells of the plate (each control must be assayed in duplicate). Add 100 µl of each sample to be tested to the remaining wells. We recommend running samples also in duplicate. Seal the plate and **incubate for 1 hour at 37°C**.
3. Throw away the liquid contained in the plate on a recipient containing NaOH, and wash 4 times following the procedure previously described.
4. Add 100 µl of conjugate I (biotin conjugate), diluted as specified before, to each well. Seal the plate and **incubate for 1 hour at 37°C**.
5. Wash 4 times following the procedure previously described.
6. Add 100 µl of conjugate II (Streptavidin-peroxidase) prepared following previous instructions to each well. Seal the plate and **incubate 30 minutes at room temperature (25°C)**.
7. Wash 5 times more following the procedure previously described.
8. Add 100 µl of the substrate to each well. Keep the plate for **10 min at room temperature**.
9. Add 100 µl of the stop solution to each well. We recommend adding this reagent following the same order in which the substrate was added.
10. Read the OD of each well at 405 nm within the following 5 min after the addition of the stop solution.

## VIII. READING AND RESULT INTERPRETATION

The reading must be made with a spectrophotometer at **405 nm**.

### **Validation of the Test:**

To validate the assay, the ratio OD (Positive Control) / OD (Negative Control) must be higher than 5.

### **Cut-off Calculation:**

The Cut-off value must be calculated as the 20% of the OD value obtained with the positive control assayed:

$$\text{Cut-off} = 0,2 \times \text{OD (Positive Control)}$$

### **Interpretation of the Results:**

When running the test in duplicate, the OD of the sample will be calculated as the arithmetic mean of the OD values in both wells.

- ⇒ All samples with OD values higher than the cut off value, must be considered as positive
- ⇒ All samples with OD values lower than the cut off value, must be considered as negative

Desarrollado y fabricado en España por:  
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tif: +34 91.368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in by: