INGEZIM PPV COMPAC

Prod Ref: 11.PPV.K3

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I. TECHNICAL BASIS

The Kit has been designed to detect, in a very easy way, specific antibodies to PPV in porcine sera samples. This Kit is based on a blocking enzymatic immunoassay. We make a brief description of the technique below:

The recombinant PPV antigen is fixed in a solid support (polystyrene plate). The serum sample is added to each well. After incubation, we add a specific monoclonal antibody (peroxidase conjugated) specific to VP2 protein of the PPV. If the serum sample does not contain antibodies against the virus, the conjugate will bind to the plate, whereas if it contains specific antibodies

to PPV, these antibodies will bind the antigen on the plate blocking the binding of the conjugate, that's why no binding of the conjugate will be appreciated.

After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled conjugate by adding the substrate that in presence of the peroxidase will develop a colorimetric reaction.

Due to the high specificity of the monoclonal used as conjugate, the assay has a high specificity and sensitivity index.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or Smoking where specimens or Kit reagents are being handle
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit use a fresh substrate preparation.
10. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.
11. Substrate and substrate buffer must be handled with care.

III. STORAGE OF COMPONENTS

Store all plates and reagents between +2°C and +8°C.

Stop solution may precipitate at this temperature. When it happens, warm at 37°C before use.

Once opened, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the 25x concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water. When ready this solution remains stable between +2°C and +8°C.

- **Preparation of the conjugate: to make immediately before use:**

Dilute the needed quantity of conjugate provided in the Kit 1/100 with diluent:

- ⇒ The necessary and sufficient quantity of conjugate for a complete plate is 60 µl of Conjugate with 6 ml of diluent.
- ⇒ The necessary and sufficient quantity of Conjugate for an eight wells strip is 5 µl of conjugate with 500 µl of diluent.

Shake very well the solution before use.

Prepare only the quantity needed for each time because the remainder volume has to be rejected.

- **Control Sera:**

Control sera must be handle as samples by adding 20 µl of serum and 80 µl of diluent.

We recommend making aliquots and storing them at -20°C.

- **Substrate**

Ready to use.

VI. TEST PROCEDURE

1. All reagents (except conjugate) must be allowed to come to room temperature before use.
2. Add 80 µl of diluent to all the wells to be used (including the ones for the positive and negative control serum).
3. Add 20 µl of Positive control serum, Negative control serum and samples to each well of the plate. For high accuracy it is recommended to run samples and controls in duplicate. Seal the plate and **incubate for 30 min at 37°C**.
4. Without removing the sera, add 50µl of conjugate, prepared as described above, to each well. Shake carefully the plate for a better homogenisation of the contents in the wells and **incubate 30 min at 37°C**.
5. Wash 4 times following the described procedure.
6. Add 100 µl of substrate to each well. Keep the plate for **10 min at room temperature**.
7. Add 100 µl of stop solution to each well. It is recommended to add this reagent in the same order than was used for adding the substrate.
8. Read the OD of each well with a spectrophotometer at 405 nm within 5 min after the addition of stop solution.

VII. READING AND RESULTS INTERPRETATION

The reading must be done with a spectrophotometer at **405 nm**.

A) Validation of the Test:

To validate the assay, the ratio OD (Negative Control) / OD (Positive control) must be higher than 4

B) Interpretation of the Results:

For each sample, the blocking % must be calculated by:

$$\frac{\text{OD (Neg. Control)} - \text{OD (sample)}}{\text{OD (Neg. Control)} - \text{OD (Pos. Control)}} \times 100$$

When you are running duplicate samples, OD values will be calculate for each sample as the arithmetic mean of both values.

- All samples with blocking % higher than 30%, must be considered as positives
- All samples with blocking % lower than 25%, must be considered as negatives
- Samples with blocking % between both values must be considered as doubtful. For these samples a new assay is recommended after 3 weeks.

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in by: