INGEZIM ROTAVIRUS PORCINO

Prod Ref: 11.RTP.K1

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I. TECHNICAL BASIS

The kit is based on an indirect enzymatic immunoassay (**Indirect ELISA**). We make a brief description of the technique below:

We fix the antigen on a solid support (Polystyrene plate). When a sample serum contains specific antibodies against the virus, they will bind to the antigen absorbed on plate. After washing to eliminate all non-fixed material from the sera sample, we can detect

the presence of porcine immunoglobulins using an specific peroxidase conjugate. After addition of the substrate, a colorimetric reaction will appear which could be measured by a spectrophotometer.

In this way, the presence of colour means the presence of antibodies against the virus in the porcine sera and the absence of colour, the absence of specific antibodies.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking or smoking where specimens or kit reagents are being handled
7. Do not pipette by mouth. Use a new tip for each serum sample.
8. For each utilisation of the kit, positive and negative control serum must be tested in a systematic way.
9. Substrate is very sensitive to light and contamination so never pipette directly from the substrate storage bottle. The necessary amounts for each assay must be always pour into a separate container.
10. Stop solution is a strong acid. Handle with care

III. STORAGE OF COMPONENTS

Plates and reagents must be stored between +2°C and +8°C.

Once opened, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

IV. INFORMATION ABOUT THE WASHING STEPS

Washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

- **Titration:**

We recommend to assay 8 two fold dilutions for each sample (from 1/200 to 1/25600). For your facilities this step could be done directly on the plate.

- **Screening:**

Sera samples must be used at 1/200 dilution in diluent to be tested (5 µl of serum in 1 ml of diluent).

VI. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionised water. When ready this solution remains stable stored at 4°C.

- **Preparation of control sera (+) and (-):**

Make a dilution 1/200 (as sample sera) in diluent (5 µl of serum in 1 ml of diluent). If you need to do a partial use of the kit. We recommend to make aliquots of the control sera after the reconstitution and to store them at -20°C till use.

- **Preparation of Diluent:**

Dilute one volume of the concentrated diluent provided in the kit with 4 volumes of distilled or deionised water. Store this solution at 4°C.

- **Preparation of the conjugate: to make immediately before use.**

Dilute 1/100 with serum diluent 1x prepared as it's said in point "Preparation of Diluent":

- The necessary and sufficient quantity of conjugate for a complete plate is 110 µl of conjugate in 11 ml of diluent.
- The necessary and sufficient quantity of conjugate for an 8 wells strip is 10 µl of conjugate in 1 ml of diluent

VII. TEST PROCEDURE:

1. All reagents (except conjugate) must be allowed to come to room temperature before use.
2. **SCREENING:**

Add 100µl of controls and samples diluted 1/200. Use two wells for positive control, two for negative control and two for each sample to be tested.

TITRATION:

Add 100 µl of each dilution of the samples (diluted as it is described in the previous instructions) using one row per sample. Add 100µl of positive control and 100µl of negative control at dilution 1/200. Use two wells for positive control and two for negative control.

Seal the plate and incubate for **30 min at 37°C**.
3. Wash 3 times following the procedure previously described.
4. Add 100 µl of conjugate, prepared following the previous instructions, to each well. Seal the plate and incubate **for 30 min at 37°C**.
5. Wash 4 times following the described procedure.
6. Add 100 µl of substrate. Keep the plate for **10 min at room temperature** in darkness.
7. Add 100 µl of stop solution to each well.
8. Read the absorbances of each well with a spectrophotometer at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND INTERPRETATION OF THE RESULTS:

The reading must be done with a spectrophotometer at 450 nm.

A) Validation of the plate:

For validation, OD values of the control sera must fulfil these requirements:

- OD of Positive control at 1/200 > 0.7
- OD of Negative control at 1/200 < 0.4

B) Results Interpretation:

The Cut Off +/- = 0.4

Samples with OD values higher than the cut off must be considered positive.

In case of **SCREENING**, the OD value of each sample will be the average between the OD values obtained in the two wells used for this serum.

In case of **TITRATION**, the titre of the sample will be the last dilution at which the OD value is higher than the cut off.

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: