## INGEZIM PRRS DR

Prod Ref: 11.PRS.K0

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## I. TECHNICAL BASIS

---

The kit has been designed to detect antibodies specific for PRRS virus, being able to detect a very low titters of antibodies in sera of infected and vaccinated animals. **The kit has been developed using the European PRRSV N protein as antigen. Nevertheless, due to the homology existing between both kinds of strains, the kit must not be used to differentiate between American and European strains. There are studies showing detection of some sera from animals infected with American strains.**

The INGEZIM PRRS DR kit is based on a novel immunoenzymatic assay called double recognition ELISA which is described below:

Plates are coated with N protein of PRRSV.

When N protein conjugated with peroxidase is added and positive sera samples contain N protein specific antibodies, they will catch simultaneously the labelled N protein and the protein N fixed to the well. Presence or absence of labelled N protein will be detected by addition of substrate (TMB) which, in presence of the peroxidase, will develop a colorimetric reaction.

This is a very sensitivity assay and it is recommended for early detection after infection/vaccination of antibodies specific of PRRSV and for negative herds certifications and testing new animals before introduced. It is not recommended for survey of vaccinated animals because all OD values will be very high.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

---

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
10. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
11. Substrate must be handle with care, it is very sensible to light and contamination.

## III. STORAGE OF COMPONENTS

---

All reagents and plates must be stored between +2°C and +8°C. **Once opened**, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

## IV. INFORMATION ABOUT THE WASHING STEPS

---

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF REAGENTS

### • **Washing solution:**

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionised water. Once prepared, this solution remains stable between +2°C and +8°C.

### • **Positive and Negative controls:**

Controls must be tested as samples (20µl / well).

### • **Preparation of the conjugate:**

Conjugate is ready to use and must be used adding 80µl/well

## VI. PREPARATION OF SAMPLES:

Sera samples must be assayed diluted 1/5 in conjugate. This dilution could be made directly on the assay plate by adding 80 µl of conjugate and 20 µl of sample to each well.

## VII. TEST PROCEDURE

- All reagents must be allowed to come to room temperature before use.
- Add in first place 80 µl of conjugate to each well and then add 20 µl of each sample to be tested.
- Add 20 µl of positive control and 20 µl of negative control to wells selected for this purpose.
- We recommend running controls and samples in duplicate. **Seal the plate and incubate 1 hour at room temperature.**
- Wash 3 times following the described procedure.
- IMPORTANT: Wash another 3 times with distilled or deionised water.**
- Add 100 µl of substrate, to each well. Keep the plate for **10 min at room temperature.**
- Add 100 µl of stop solution to each well.
- Read the OD of each well at 450 nm within 5 min after the addition of stop solution.

## VIII. READING AND RESULT INTERPRETATION

Determine the mean absorbance of positive and negative controls and samples, in case of assayed the last ones in duplicate.

### 1. VALIDATION OF THE RESULTS

OD of Positive Control > 1,4  
OD of Negative Control < 0,25

### 2. CUT OFF CALCULATION

Cut Off = OD of positive control x 0.175

### 3. RESULTS INTERPRETATION

Samples will be considered **POSITIVE** (there are antibodies specific of PRRSV), if the OD value at 450nm is higher than the cut off.

Samples will be considered **NEGATIVE** (there are no antibodies specific for PRRSV in the sample) if the OD value at 450nm is equal or lower than the positive cut off.

Because this assay is highly sensitivity it is recommended with weak positive animals to take another sample 5-7 days after. If the animal is a recent infected, this samples will be a high positive.

**ATTENTION: The kit has been developed using the European PRRSV N protein as antigen. Nevertheless, due to the homology existing between both kinds of strains, the kit must not be used to differentiate between American and European strains.**

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: ingenasa@ingenasa.com  
www.ingenasa.es

Distributed in by: