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
Enzymatic immunoassay for detection and quantification of antibodies specific for porcine respiratory and reproductive syndrome virus (European and American strains) in swine serum

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OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT:

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Distilled or deionised water.
Micropipets from 5 to 200 µl.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 ml
ELISA Reader (405 nm filter)

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I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). A brief description of the technique is given below:

The antigen is coated on polystyrene plates. When a sample serum contains specific antibodies against the virus, they will bind to the antigen adsorbed on the plate. After washing to eliminate all non-fixed material the presence of swine antibodies can be detected using a specific peroxidase conjugate. After addition of the substrate a colorimetric reaction will appear that can be measured by a spectrophotometer. The presence of colour indicates the presence of antibodies against the virus in the swine sera, and the absence of colours, the

absence of specific antibodies. The aim of this kit is to provide the users with a reliable and automatizable diagnostic technique for this disease, for which the immunoperoxidase monolayer assay and immunofluorescence had to be used till now.

Please note that our kit uses purified antigen obtained by the expression of the ORF7 from American and European strains of PRRS virus in E.Coli. This method warrants the total absence of infectivity in the kit.

The comparison assay made between this technique and the immunoperoxidase and immunofluorescence assays gives very satisfactory correlation levels.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions for use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiry dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Substrate and substrate buffer must be handled with care.
10. For each utilisation of the Kit use a fresh substrate preparation.
11. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.

III. STORAGE OF COMPONENTS

Store all plates and reagents between +2°C and +8°C.

WARNING! Stop Solution may precipitate, if this happens, warm the solution to +37°C prior to use.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.
- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.

V. PREPARATION OF SAMPLE

Samples to be used in the assay must be porcine sera at **1/100 dilution in diluent DE13-01**

(i.e.: 5 µl of serum in 0,5 ml of serum diluent DE13-01).

VI. PREPARATION OF REAGENTS

- **Washing solution:**
Dilute one part of the concentrated washing solution provided in the kit with 24 parts of distilled or deionized water (40 ml of concentrated solution and 960ml of water). Once ready this solution remains stable at +4°C.
- **Diluents:**
Both diluents are supplied ready to use. **Do not dilute.**
- **Control sera:**
Control sera are ready to use (1/100 dilution). **Do not dilute.**
- **Preparation of the conjugate: prepare immediately before use.**
Dilute 1/100 the required quantity of conjugate provided in the kit with diluent DE01-01 (e.g. 110 µl of conjugate into 11 ml of conjugate diluent is sufficient for one 96- well plate). Shake the solution thoroughly before use. Prepare only the quantity required for each use as any surplus solution will need to be discarded.
- **Substrat:**
Ready to use.

VII. TEST PROCEDURE

1. All reagents (except conjugate) must be allowed to reach room temperature before use.
2. Add 100 µl of positive control to two wells of the plate, 100 µl of the negative control to another 2 wells and 100 µl of sera samples to the remaining wells of the plates (for the purpose of confirmation we recommend to run samples in duplicate wells). Seal the plate and incubate for 1 h at 37°C.
3. Wash 4 times following the procedure previously described.
4. Add 100 µl of conjugate, prepared following the previous instructions to each well. Seal the plate and incubate for 45 minutes at room temperature (20-25° C)
5. Wash 5 times following the described procedure.
6. Add 100 µl of substrate to each well. Keep the plate at room temperature for 15 minutes.
7. Add 100 µl of stop solution to each well.
8. Read the OD of each well at 405 nm within 5 min after the addition of stop solution.

VIII. READING AND INTERPRETATION OF THE RESULTS:

Validation of the test:

The test is considered valid when positive control OD is higher than 1.500 and negative control OD is lower than 0.250.

Interpretation of the results:

1. CUTT OFF CALCULATION:

The cut off must be calculated as follows:

$0.15 \times \text{OD of positive control.}$

Samples with OD higher than cut off value should be considered positive to PRRS antibodies.

Samples with OD lower than cut off value should be considered negative to PRRS antibodies.

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The titre of positive samples will be determined as follows:

- a) Calculate S/P

$$S/P = \frac{\text{Sample OD}}{\text{Positive Control OD}}$$

b) Titre is calculated as follows

$$\text{Titre} = 100 (e^{4x}) - 80$$

Where,

e is the natural logarithm base and **x** is the S/P of the sample.

EXAMPLE:

Serum	OD
Control +	2.343
Control -	0.182

Sample	OD
Sample "1"	0.910

1. **Validation:**

Positive control OD value = 2.343 (> 1.5)

Negative control OD value = 0.182 (< 0.25)

2. **Cut off calculation**

$$0.15 \times 2.343 = 0.351$$

$$\text{Sample OD} = 0.910 \rightarrow \rightarrow \text{positive}$$

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a) $S/P (x)$

$$x = 0.910 / 2.343 = 0.3884$$

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$$\text{Titre} = 100 (e^{4 \times 0.3884}) - 80 = 393$$

All calculations can be done using Microsoft Excel software

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in by: