## INGEZIM M.HYO COMPAC

Prod Ref: 11.MHY.K3

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## I. TECHNICAL BASIS

---

The kit have been designed to detect antibodies specific for Mycoplasma hyopneumoniae (M.hyo) in swine serum.

The INGEZIM MHYO kit is based on the blocking immunoenzymatic assay described below:

Plates are coated with M.hyo antigen. After adding the sample to the well, if it contains specific antibodies against M.hyo, they will bind to the antigen absorbed on plate while if the sample does not contain specific antibodies they will not bind to the antigen.

If we add a specific monoclonal antibody against the M.hyo antigen coated to the plate (conjugated with peroxidase), it will compete with the antibodies of the serum. If sera samples contain specific antibodies, they will not permit the binding of the labelled MAb to the antigen whereas if it does not contain specific antibodies the MAb will bind to the antigen on the plate. After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled MAb by adding the substrate (TMB) that in presence of the peroxidase will develop a colorimetric reaction.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

---

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
10. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
11. Substrate must be handle with care, it is very sensible to light and contamination.

## III. STORAGE OF COMPONENTS

---

All reagents and plates must be stored between +2°C and +8°C.

## IV. INFORMATION ABOUT THE WASHING STEPS

---

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF REAGENTS

### • Washing solution:

Dilute one part of the 25x concentrated washing solution provided in the Kit with 24 parts of distilled or deionised water (i.e. 40 ml of concentrate and 960 ml of water). Once prepared, this solution remains stable between +2°C and +8°C.

### • Positive and Negative controls:

Controls are supplied ready to use (Add 100µl/well)

### • Preparation of the conjugate

**CONJUGATE IS SUPPLIED READY TO USE: DO NOT DILUTE**

## VI. PREPARATION OF SAMPLES:

Sera samples must be assayed at a 1/2 dilution in diluent. This dilution could be made directly on the assay plate by adding 50 µl of diluent and 50 µl of sample to each well.

## VII. TEST PROCEDURE

- All reagents must be allowed to come to room temperature before use.
- Add 100µl of positive control to two wells and 100 µl of negative control to another two wells. Add 50 µl of diluent to remainder wells and then, add 50 µl of each sample to be assayed. We recommend running controls in duplicate wells. Seal the plate and incubate 1 hour at room temperature (18-25°C).
- Wash 4 times following the described procedure.
- Add 100 µl of conjugate, to each well. Seal the plate and incubate for 30 min at room temperature (18-25°C).
- Wash 6 times following the described procedure.
- Add 100 µl of substrate, to each well. Keep the plate for 10 min at room temperature.
- Add 100 µl of stop solution to each well.
- Read the OD of each well at 450 nm within 5 min after the addition of stop solution.

## VIII. READING AND RESULT INTERPRETATION

In case that the samples had been run in duplicate, it has to be considered the mean of both OD values. In the same way, the mean of the values obtained in the two well of positive and the two well of negative control has to be made.

### 1. VALIDATION OF THE RESULTS

OD Positive Control / OD Negative Control < 0.25  
OD Negative Control > 0.75

### 2. CUT OFF CALCULATION

Cut Off (-) = 0.45 x Negative Control

Cut Off (+) = 0.40 x Negative Control

### 3. BLOCKING % CALCULATION

$$\text{Blocking \% of sample} = 100 - \frac{\text{Sample OD} \times 100}{\text{Negative control OD}}$$

### 4. RESULTS INTERPRETATION

Samples will be considered **POSITIVE** (there are antibodies specific for M.hyo), when de OD at 450nm was equal or lower than the positive cut off (40 % of negative control).

Samples will be considered **NEGATIVE** (there are no antibodies specific for M.hyo in the sample) when the OD value at 450nm was equal or higher than the negative cut off (45% of negative control).

Samples with percentages of binding between both cut off have to be considered **DOUBTFUL**. In these cases a new sample of the animal is recommended to be analysed in two weeks.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tif: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: ingenasa@ingenasa.com  
www.ingenasa.es

Distributed in \_\_\_\_\_ by: