INGEZIM APP mix

Prod Ref: 11.APP.K1

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I. TECHNICAL BASIS

This is an immunoenzymatic assay for the detection of antibodies against APP parasuis in porcine serum. The porcine serum samples and the controls are diluted and incubated in wells coated with APP antigens. The antibodies (Abs) specific to APP present in positive serum samples will bind to the Ags in the wells. After several washes to eliminate unbound substances, a conjugate (an Ab coupled to an enzyme) targeted at the porcine Abs is added. After

incubation, the excess of this conjugate is eliminated by a second wash and its attachment is revealed with a chromogenous substrate. Following this incubation, the enzyme, if present, reacts with the substrate and a green color develops. The reaction is stopped and the optical densities are read. The intensity of the colour allows the determination of the type of sample tested.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.
10. Substrate and stop solution must be handled with care.

III. STORAGE OF COMPONENTS

All reagents and components must be stored between +2°C and 8°C. WARNING: Stop Solution may precipitate at this temperature, if it happens, warm up at +37°C before use.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLE

It is strongly recommended to test the samples and controls as duplicate.

Dilute porcine serum samples at 1/200 in 1X serum diluent (e.g., 4 µl sample in 796 µl 1x diluent). Make sure you use a new tip for each sample. Also make sure each dilution is properly mixed before being distributed into the wells.

VI. PREPARATION OF REAGENTS

• **Washing solution:**

Dilute one part of the concentrate washing solution provided in the kit with 24 parts of distilled or deionized water (i.e. 40 ml of concentrate with 960 ml of water).

• **Preparation of Diluent:**

Supplied diluent is 10x concentrated. Before using it, a 1/10 dilution must be done with distilled or deionized water (1 vol. of concentrated supplied Diluent with 9

volume of distilled water). When ready this solution remains stable at +4°

• **Control sera:**

Controls are supplied ready to use. Do not dilute

• **Preparation of the conjugate: to make immediately before use.**

Dilute the needed quantity of conjugate provided in the kit 1/100 in diluent:

The necessary and sufficient quantity for a complete plate is 110 µl of conjugate with 10,9 ml of diluent.

Shake very well the solution before the use. Prepare only the quantity needed for each time because the remainder volume has to be rejected.

• **Preparation of the substrate solution.**

Substrate is supplied ready to use.

VII. TEST PROCEDURE

Bring all reagents to room temperature and mix well manually before use.

1. Make a schematic representation of the plate and the distribution of controls and samples.
2. Dispense 100 µL of each ready-to-use positive control into two wells:
 - For APP 1-2-9-11 plates, use controls APP 1-9-11 and APP 2
 - For APP 3-6-8-15 plates, use control APP 3-6-8-15
 - For APP 4-5-7 plates, use controls APP 4-7 and APP 5
3. Dispense 100 µL ready-to-use negative control into two wells of each of the plates.
4. Dispense 100 µL diluted samples (see section B) into two wells of each of the plates APP 1-2-9-11, APP 3-6-8-15, APP 4-5-7.
5. Incubate at $23 \pm 2^\circ\text{C}$ for 30 minutes.
6. Wash each well 5 times with 300 µL 1X wash solution (see section A). Throw away all liquid

7. Dispense 100 µL diluted conjugate (see section C) into each well.
8. Incubate at $23 \pm 2^\circ\text{C}$ for 30 minutes.
9. Repeat step 6.
10. Dispense 100 µL ready-to-use substrate into each well.
11. Incubate, away from light, at $23 \pm 2^\circ\text{C}$ for 20 minutes.
12. Dispense 100 µL ready-to-use stop solution into each well.
13. Measure optical densities (OD) at 405 nm. If the microplate reader is equipped with a reference filter, set it at 490 nm. The reading should be done no later than 15 minutes after the addition of the stop solution.
14. Calculate the results.

VIII. READING AND RESULT INTERPRETATION

1. Calculation of the Ratio

APP 1-2-9-11

Calculate the mean (M) of ODs obtained for the APP 1-9-11 positive control (M OD APP 1-9-11), the APP 2 positive control (M ODAPP 2), the negative control and each sample (M ODS). Then, calculate the mean of the means of the positive controls (M ODAPP 1-9-11-2). To obtain ratio, divide each sample's mean OD by the mean of the means of the positive controls (M OD APP 1-9-11-2).

Example

- OD APP 1-9-11 : 1.250 et 1.200 ; M ODAPP 1-9-11 = $(1.250 + 1.200)/2 = 1.225$
- OD APP 2 : 1.100 et 1.200 ; M OD APP 2 = $(1.100 + 1.200)/2 = 1.150$
- o M OD APP 1-9-11-2 = $(1.225 + 1.150)/2 = 1.187$
- OD S = 0.800 et 0.850 ; M OD S = $(0.800 + 0.850)/2 = 0.825$
- Ratio S = $0.825/1.187 = 0.695$

APP 4-5-7

Proceed as for plates APP 1-2-9-11

APP 3-6-8-15

Calculate the mean (M) of ODs obtained for the APP 3-6-8-15 positive control (M OD 3-6-8-15), the negative control and each sample (M ODS). To obtain ratio, divide each sample's mean OD by the mean value of the positive control (M OD APP 3-6-8-15).

2. Validity Criteria

The following criteria must be met in order to validate the analysis:

Negative control ratio must be less than 0.15.

Mean of each of the different positive controls ODs must be greater than 0.70.

3. INTERPRETATION

PLATES APP 1-2-9-11

- RATIO < 0.33: Negative sample
- RATIO \geq 0.5: Positive sample
- $0.33 \leq$ ratio < 0.5: Suspicious samples

PLATES APP 3-6-8-15

- RATIO < 0.30: Negative sample
- RATIO \geq 0.5: Positive sample
- $0.30 \leq$ ratio < 0.5: Suspicious samples

PLATES APP 4-5-7

- RATIO < 0.35: Negative sample
- RATIO \geq 0.5: Positive sample
- $0.35 \leq$ ratio < 0.5: Suspicious samples

In case of suspect or positive result, it is recommended to test the samples against the relevant serotypes using kits specific for these serotypes (ie. Swinecheck® APP 1-9-11, APP 2, APP 3-6-8, APP 4-7, APP 5).

Desarrollado y producido por:
Developed and manufactured by:

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INMUNOLOGIA Y GENETICA APLICADA, S.A.