INGEZIM ANEMIA EQUINA

Prod Ref: 14.AIE.K0

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I. TECHNICAL BASIS

The kit has been designed to detect antibodies specific for Equine Infectious Anemia Virus (EIAV) in horses being able to detect a very low titters of antibodies in sera of infected animals.

The INGEZIM EIA DR kit is based on a double recognition ELISA which is described below:

Plates are coated with recombinant P26 protein of EIAV. After adding the sample to the well, if it contains EIAV specific antibodies, they will bind to

the antigen coating the plate while if the sample does not contain specific antibodies they will not bind to the antigen.

When P26 protein conjugated with peroxidase is added and sera samples contain P26 specific antibodies, they will catch the labelled P26 while they are fixed to the antigen coated in the well. Presence or absence of labelled P26 will be detected by addition of substrate (TMB) which, in presence of the peroxidase, will develop a colorimetric reaction.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
10. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
11. Substrate must be handled with care, it is very sensible to light and contamination.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C. **Once opened**, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.

Dispense a volume of 300 µl of washing solution on each well.

Shake delicately the plate, avoiding the contamination between wells.

Turn over the plate brusquely to empty the wells.

Repeat the process as much times as is indicated on the instructions of the Kit.

Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.

After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF REAGENTS

• Washing solution:

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionised water (i.e. 40 ml of concentrate and 960 ml of water). Once prepared, this solution remains stable between +2°C and +8°C.

▪ Positive and Negative controls:

Controls must be tested as the samples (adding 50µl diluent + 50 µl control)

▪ Preparation of the conjugate:

Conjugate is ready to use and must be used adding 100µl/well

VI. PREPARATION OF SAMPLES:

Sera samples must be assayed diluted 1/2in diluent. This dilution could be made directly on the assay plate by adding 50µl of diluent and 50 µl of sample to each well. **ATTENTION: HAEMOLYZED OR CONTAMINATED**

SAMPLES MAY GIVE FALSE RESULTS

VII. TEST PROCEDURE

- All reagents must be allowed to come to room temperature before use.
- Add 50 µl of diluent to each well. Add 50 µl of each sample to be assayed and then, 50 µl of positive and negative control. We recommend running controls in duplicate wells. For confirmatory purposes is recommended but not mandatory to run samples in duplicate wells Seal the plate and incubate 1 hour at 37°C.
- Wash 3 times following the described procedure.
- Add 100 µl of conjugate ready to use to each well. Seal the plate and incubate for 1 hour at 37°C.
- Wash 3 times following the described procedure.
- Add 100 µl of substrate, to each well. Keep the plate for 15 min at room temperature.
- Add 100 µl of stop solution to each well.
- Read the OD of each well at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

If you run samples and controls in duplicate wells calculate the mean OD of both wells.

1. VALIDATION OF THE TEST

OD of Positive Control > 0.8

OD Negative Control/ OD Positive Control < 0.15

2. RESULTS INTERPRETATION

Calculate the relation OD of sample/OD of positive control (S/P).

Samples will be considered **POSITIVE** (there are antibodies specific of EIAV), if the S/P value is higher or equal than 0.3.

Samples will be considered **NEGATIVE** (there are no antibodies specific for EIAV in the sample) if the S/P value is lower than 0.3.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by: