## INGEZIM AIE-CROM

Prod Ref: 14.AIE.K41

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## I. TECHNICAL BASIS

---

The test is based on the technique of immunochromatography, a migration technique that uses purified RP26 protein of Equine Infectious Anaemia (EIA), and is capable of detecting specific antibodies against EIA in horse serum samples.

The diagnostic device consists of two windows:

- Sample window: it contains RP26 protein of EIA and one protein control, coated to coloured latex particles.
- Result window: it contains a test line (T) formed by the RP26 protein and a control line (C) formed by a MAb specific for the protein control.

If the sample contains antibodies against the RP26, they will react with the red particles conjugated to RP26 protein. The latex-protein-antibody complex migrates through the membrane and reacts again with the RP26 protein adsorbed on the test line resulting in the appearance of a red/pink. The appearance of a blue line in the control area (C) indicates that the chromatography has been correctly performed.

## II. STORAGE:

---

All the reagents must be kept between 4 and 25 °C.

## III. PRECAUTIONS AND WARNINGS FOR USERS:

---

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handle

## IV. TEST PROCEDURE

---

Reagents and strips must be at room temperature when used for testing.

1. Open the aluminium bag just in the moment to realise the assay. Remove the diagnostic device from the bag and place it on a flat surface.
2. Dispense 10 µl of undiluted serum at the round window of the diagnostic device.
3. Slowly, add 3-4 drops of running buffer on the same round window.
4. Read the results at 10 minutes.


## V. INTERPRETATION OF RESULTS

---


In order to be sure that the chromatographic assay has run properly, a blue line must appear in control area (C).

- **POSITIVE:** Two lines appear in the window reading: a blue line (control) and a pink/red (test line).
- **NEGATIVE:** only a blue line appears (control).
- **NO VALID TEST:** No blue line appears (control). In this case, the test can not be considered valid and should be repeated with a new strip.

No results must be considered after 10 minutes.

1. Añadir 10  $\mu$ l de suero en la ventana circular  
Add 10  $\mu$ l of serum sample on the round window

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INMUNOLOGIA Y GENETICA APLICADA, S.A.