INGEZIM AHSV COMPAC PLUS

Prod Ref: 14.AHS.K3

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by:

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Blocking immunoenzymatic assay, for
detection of specific antibodies to African
Horsesickness Virus in equine sera.

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Distilled or deionised water.
Micropipettes from 5 to 200 µl.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 ml
ELISA Reader (450 nm filter)

VI. PREPARATION OF REAGENTS

- **Washing solution:**
Dilute one part of the concentrate washing solution provided in the Kit in 24 parts of distilled or deionized water (40 ml of concentrate solution and 960 ml of water). Once this solution is ready, it remains stable between +2°C and +8°C.
- **Preparation of conjugate:**
Conjugate is ready to use and does not need any preparation.
- **Preparation of substrate:**
Ready to use
- **Preparation of Controls (+) y (-):**
Controls are ready to use and do not need any preparation. Dispense 100 µl of each one.

VII. TEST PROCEDURE

1. All reagents must be brought to room temperature before use.
2. Dispense 100 ul of diluted samples into appropriate wells (dilution 1/5). We recommend running the samples in duplicate. Dispense 100 ul of positive control into two wells and 100 ul of negative control into two wells. Cover the plate and **incubate for 1 hour at 37°C.**
3. Wash 5 times as described before
4. Add 100 ul/well of conjugate. **Incubate for 30 minutes at 37°C.**
5. Wash 5 times as described before.
6. Dispense 100 ul/well of substrate using a multichannel pipette. **Incubate for 10 minutes at room temperature.**
7. Dispense 100 ul/well of stop solution. Take care to avoid dispensing of bubbles.
8. Read at 405nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **405 nm.**

A. Validation criteria:

- The OD of the Positive control must be lower than 0.2
- The OD of the Negative control must be higher than 1.0

B. Interpretation:

Determine the Blocking Percentage (BP) of each sample applying the following formula:

$$BP = \frac{\text{Abs}(\text{Control -}) - \text{Abs}(\text{sample})}{\text{Abs}(\text{Control -}) - \text{Abs}(\text{Control +})} \times 100$$

D. Example

Average Abs Positive Control = 0,110
Average Abs Negative Control = 1,859
Abs sample = 0,865

$$BP = \frac{1,859 - 0,865}{1,859 - 0,110} \times 100 = 56,8$$

Sample is positive.

I. TECHNICAL BASIS

This kit is based on a blocking enzymatic immunoassay (Blocking Elisa). A brief description of the technique is shown below:

The antigen is fixed in a solid support (polystyrene plate). After incubation with sera samples, an AHSV specific monoclonal antibody (Mab peroxidase conjugated) is added. If the sample contains antibodies specific of the virus, they will not allow the binding of the labelled Mab to the antigen whereas if it does not contain specific antibodies, the Mab will bind to the antigen coating the plate.

After washing the plate to eliminate all non-fixed material, presence or absence of labelled Mab can be detected by adding the substrate which, in presence of the peroxidase, will develop a colorimetric reaction.

If there has been colour development, the conjugate has bound to the antigen, being the sample negative. On the other hand, if there are antibodies specific of AHSV in the sample, they will block the binding of the conjugate and there will not be colour development. The antigen used in this kit is VP7 recombinant protein from the AHSV (serotype 4) obtained using the baculovirus expression system. This confers many advantages to the assay. From the safety point of view, it is clear the total absence of infectivity. From the accuracy point of view, VP7 is one of the mayor proteins from AHSV, as well as the most antigenic one and the most conserved within the 9 different serotypes. It has not been found any infected or vaccinated animal without antibodies to this protein.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the kit's reagents.
5. Do not use components after expiration dates and do not mix components from different lots.
6. Do not eat, drink, or smoke where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each use of the kit, control positive and negative must be tested in a systematic way.
10. Substrate must be handled with care, it is very sensible to light and contamination.
11. Stop solution can precipitate at refrigeration. Be sure there is not precipitation before use

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

The stop solution precipitates at this temperature. Incubate at 37°C before use, to eliminate the precipitation.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

Sera samples must be diluted 1/5 in the provided diluent. This step can be made directly in the wells of the plate by adding 80 µl of diluent and 20 µl of serum sample and shaking carefully the plate to homogenize it.

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