INGEZIM PEA DAS

Prod Ref: 14.PEA.K2

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INGEZIM PEA DAS 14.PEA.K2

I. TECHNICAL BASIS

This kit is based on a double antibody sandwich enzymatic immunoassay (DAS or capture Elisa). We make a brief description of the technique below:

The plate is coated with a specific antibody to vp7 of the AHSV. In this way when a sample is added into the wells, if it contains antigen, the antibodies

will catch it. After incubation and washing a biotin-conjugated antibody to vp7, is added. If there is antigen captured into the well, the biotin conjugate will bind to it. After incubation and washing you add a second conjugate (Streptavidin peroxidase conjugated) and the peroxidase action will be detected after adding the adequate substrate.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the kit's reagents.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
10. Substrate must be handled with care, it is very sensitive to light and contamination.
11. Stop solution can precipitate at refrigeration. Be sure there is not precipitation before use.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

The stop solution is precipitated at this temperature. Incubate at 37°C before using, to eliminate the precipitation.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

Any kind of biological samples could be assayed in this kit. Nevertheless, for diagnosis purposes, we recommend to use the **spleen**

To prepare the sample, 1/2 dilution must be done (w/v) using the diluent supplied (i.e. 1g of

sample in 2 ml of diluent). For optimal homogenisation of the sample use a stomacher or any other mechanic device. After that, centrifuge at 1000 rpm for 10 minutes and take the supernatant for testing.

VI. PREPARATION OF REAGENTS

- **Washing solution:**
Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water (40 ml of concentrate solution and 960 ml of water). Once this solution ready, it remains stable between +2°C and +8°C.
- **Controls:**
Controls are ready to use. (Add 200 µl / well)
- **Preparation of the conjugates: to make immediately before use.**

Dilute the needed quantity of conjugates provided in the kit 1/100 with diluent (e.g. 20 µl of conjugate into 2 ml of diluent is enough quantity for one strip of 8 wells). Shake the solution before use.

ATTENTION: Prepare only the quantity needed for each time because the remainder volume has to be rejected.

- **Preparation of the substrate:**
Ready to use.

VII. TEST PROCEDURE

1. All reagents (except conjugates) must be allowed to come to room temperature before use.
2. Add 200 µl of positive control to one well (e.g. A1), and 200 µl of the negative control to other well (e.g. A2). Add 200 µl of samples to test on each remainder wells. We recommend the use of two wells per sample. Seal the plate and incubate **for 1 h at 37°C**.
3. **Empty the wells into a receptacle containing 0.1M NaOH** and wash four times as it is specified at point IV.
4. Add 200 µl of biotin conjugate (prepared following previous instructions) to each well. Seal the plate and **incubate for 1 hour at 37°C**.
5. Wash 5 times following the described procedure.
6. Add 200 µl of peroxidase conjugate (prepared following previous instructions) to each well. Seal the plate and **incubate for 45 minutes at room temperature (20-25°C)**.
7. Wash 5 times following the described procedure.
8. Add 200 µl of substrate solution, prepared following previous instruction, to each well. **Keep the plate for 10 minutes at room temperature.**
9. Add 50 ul/well of stop solution
10. Read a 405 nm within 5 min after the addiction of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **405 nm**.

A. Validation criteria:

The test will be considered valid when:

- The OD of the Positive control is higher than 2.0

B. Results Interpretation:

If sera are assayed in duplicates, OD of the sample will be calculated as the arithmetic mean of OD

values in both wells. The corrected OD value should be calculated for each sample:

Corrected OD = Sample OD – Negative Control OD

- Samples with a corrected OD higher than 0.200 must be considered as **POSITIVE** to AHSV.
- Samples with corrected OD lower than 0,150 must be considered as **NEGATIVE** to AHSV.
- Samples with OD values between both values are considered as **DOUBTFUL**. We recommend re-testing these animals applying a different technique.

IX. IMPORTANT NOTE:

This assay has been checked with spleen only for serotype 4. The rest of serotypes have been checked with cell culture supernatants.

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by: